



way

catalog & manual












Index

| | |
|----|---|
| 5 | Key |
| 6 | Way implant system |
| 8 | Packaging |
| 9 | Color code |
| 10 | Surgical protocol |
| 10 | Preparation of the implant site |
| 15 | Removal of implant |
| 16 | Implant insertion |
| 18 | Tightening of the cover screw |
| 19 | Instruments |
| 19 | Surgical organizers |
| 24 | Drills |
| 30 | Spanners and inserts |
| 36 | Surgical planning |
| 36 | Implantology equipment |
| 38 | way Mix |
| 40 | way Extra |
| 42 | Prosthetic components way Mix and way Extra |
| 52 | way Slim |
| 53 | Prosthetic components way Slim |
| 56 | way Rock |
| 58 | way Short |
| 59 | Prosthetic components way Rock e way Short |
| 65 | Digitals |
| 66 | Geadrive |
| 71 | Performa |
| 74 | Advanced |
| 81 | Warnings and sales conditions |

Index

Material

| | | | | | |
|--|------------------------|---|------------------------|---|-------------------------|
|  Inox | Stainless steel |  PEEK | Polyetheretherketone |  PTFE | Polytetrafluoroethylene |
|  Al | Aluminium |  PA | Polyamide |  PP | Polypropylene |
|  WC | Tungsten carbide |  PPSU | Polyphenylsulphone |  PU | Polyurethane |
|  CoCr | Cobalt chrome |  PMMA | Polymethylmethacrylate |  SI | Silicone |
|  EVA | Ethylene vinyl acetate |  POM | Polyoxymethylene |  Ti | Titanium |
|  Au lega | Gold alloy |  PPS | Polyphenylene sulfide | | |

Handling instruments

| | | |
|--|---|---|
|  microsam |  equator |  stepper |
|--|---|---|

Restorative components

| | |
|---|---|
|  non rotating |  rotating |
|---|---|

Abbreviations

| | | | |
|----------|----------|-----------|----------|
| H | height | L | length |
| ø | diameter | p. | platform |

The measurements shown in the catalog are expressed in mm.
The images shown are exclusively representative of the products.

way

Way is an implant system designed by Geass to answer all the professionals' needs, uniting surgical and restorative simplicity and freedom: five types of implant, specific for every kind of rehabilitation, connected by the same surgical protocol.

Way guarantees elevated levels of safety in all phases of the implant-restoration treatment, as each component is produced respecting the highest quality standards to ensure precision and functionality.



► Implant Guarantee

A guarantee for reliability that Geass is able to offer for life on a range of implants and components thanks to:



- technical reliability from thirty years of know how;
- clinical reliability from a protocol which has been applied for over 20 years.

► Synthegra surface

Synthegra is the innovative implant surface patented by Geass and applied to all the way implants, created to satisfy two requirements at the same time: to reduce the risk of peri-implant infections and to promote osseointegration.

Synthegra offers extraordinary results because it has an extremely smooth nature, able to obstacle bacterial adhesion, which at the same time acts like a rough surface, favoring osseointegration.

Synthegra is the safe and avantgarde answer to reduce the risk of peri-implantitis and to achieve the new challenge of long term osseointegration.





way Mix (p. 38)

The implant which acts on the key factors of the esthetic result: maintenance of bone levels, effective conditioning of the soft tissues, creation of an esthetically guided restoration.



way Extra (p. 40)

Designed for the rehabilitation of post-extractive sites with contemporaneous placement of the implant: extra aggressive, extra stable.



way Slim (p. 52)

The implant with 3 mm diameter, designed for reduced anatomical spaces such as the inferior incisors and the upper laterals, not treatable with traditional diameters.



way Rock (p. 56)

The transmucous implant for the management of the specifics of distal sectors, favoring restoration.



way Short (p. 58)

The implant with a length of 5 and 6.5 mm, specifically for the treatment of cases with reduced bone height, avoiding regeneration interventions.




Packaging

The packaging of the implants is characterized by:

- blister in PETG and Tyvek® to guarantee sterility;
- informative label placed on the back of the blister
- sealing sticker which, as well as guaranteeing that the packaging is not damaged, has a color code for the diameter and length of the implant.



► Information label




► Sterilization


The implants are sterilized with ionizing radiation according to the protocol validated based on current regulations. All of the other products are supplied decontaminated in non sterilizable packaging.

Color code

way|mix




way|extra




way|slim



way|rock



way|short



Prosthetic components

unique connection



small
profile



large
profile

Prosthetic components



p.3.0

Prosthetic components



connection
p.4.8



connection
p.6.5

Preparation of the implant site

The modalities and instruments for the preparation of the implant site are the same for all of the way lines and only depend on the diameter of the implant and the type of bone*.

D1

Dense cortical bone which requires further widening of the site for some implant lengths. After the standard sequence, the surgical protocol requires the use of a final drill of the same diameter of the implant but shorter, to be sunk until the definitive depth.



D2-D3

Dense/porous cortical bone, with a tight trabecular structure, or thin and porous cortical bone, with a sparse trabecular structure.



D4**

Practically non-existent cortical bone, with a sparse trabecular structure which requires the site to be underprepared. The surgical sequence requires that the last final drill is the shortest available.




*The classification is the one created by Misch (Bone character: second vital implant criterion, Dent Today 7:39-40,1998), which distinguishes four types of bone density based on the macroscopic characteristics of the cortical and trabecular bone of the edentulous portion to be treated.

** The use of way Slim and way Short implants is not foreseen in cases of D4 bone.

► Surgical sequence Ø 3



way Slim




In case of compact bone, after the standard sequence widen the perforation with the enlarger drill, sinking it up to the mark.

► Surgical sequence ø 3.4




way Mix way Extra way Rock



In case of bone D4, after the twist drill 2.1 the last step foresees to use the final drill 3.4 L. 9 mm, whatever the implant length is.


► Surgical sequence ø 3.8



In case of compact bone, after the standard sequence, widen the hole with a final drill 3.8 of inferior length (see table), sinking it up to the final depth of the site. For implants L. 5 and 6.5 mm, insert the implant for some threads; when the bone compression is too high, proceed with tapping with the implant itself, unscrewing it for a few turns and then screwing it again. Proceed with this movement until reaching the desired depth.

In case of bone D4, after the twist drill 2.5 the last step foresees the use of the final drill 3.8 L. 8 mm, whatever the implant length is

► Surgical sequence ø 4.5




In case of compact bone, after the standard sequence, widen the hole with a final drill 4.5 of inferior length (see table), sinking it up to the final depth of the site. For implants L. 5 and 6.5 mm, insert the implant for some threads; when the bone compression is too high, proceed with **tapping** with the implant itself, unscrewing it for a few turns and then screwing it again. Continue with this movement until reaching the desired depth.

In case of bone D4, after the final drill 3.8 the last step foresees the use of the final drill 4.5 L 8 mm, whatever the implant length is.

► Surgical sequence ø 5.5



way Short



osso D1 In case of compact bone insert the implant for some threads; when the bone compression is too high, proceed with **tapping** with the implant itself, unscrewing for a few turns and then screwing it again. Continue with this movement until the desired depth.



Removal of implant

► Touch&go

The touch&go functional solution is an innovative system which allows for the removal of the implant in a rapid and sure fashion without compromising its sterility. Its special ergonomics allows you to block the implant in place, facilitating coupling between the implant seat and the insert.



Before opening the implant packaging, check on the label on its back that the diameter and length measurements of the implant are suitable to the intervention. Opening of the blister must be carried out according to the clinician's own procedure to maintain sterility.



Keep touch&go in a vertical position and remove the upper part which contains the cover screw.



Press the extruding parts so that the two titanium sheets move towards each other, always keeping touch&go in a vertical position; in this way the implant is stable.

With the other hand, introduce the driver or the insert in the implant seat and match the two devices, slightly.

Reduce the pressure on touch&go and remove the implant.



Remove the cover screw contained in the upper body of the touch&go using the driver or Microesam terminal.




Implant insertion



Manual insertion




Remove the implant from the touch&go holder using the W-Start screwdriver, identifiable by the presence of the o-ring. Insert the implant for a few threads into the implant site.



Complete the insertion of the implant at crest level, using the ratchet wrench with the W-Fix insert. do not exceed the torque of 50 Ncm. Verify through the hexagons of the inserts that the orientation of the implant seat promotes the correct use of abutments.


Insertion with micromotor



Remove the implant from the touch&go holder using the W-Start driver, identifiable by the presence of the o-ring. Insert the implant, keeping below 15 rpm; do not exceed the torque of 50 Ncm.

Way Extra is indicated exclusively for rehabilitations which foresee the insertion of the implant at the same time as dental extraction (post extractive sites). In these clinical situations, it is necessary to take some other aspects into consideration:


- it is necessary to confirm the presence of 3-4 mm of native bone apically to the alveolar, essentially for the retention of the implant;
- to achieve excellent esthetic results, it is recommended that you position the implant 1 to 3 mm below the crestal margin;
- the way Extra implant makes it possible to slightly modify the insertion axes during placement, in order to achieve an excellent prosthetic orientation.





way Rock way Short


Manual insertion



Remove the implant from the touch&go holder, using the Beak insert.

Do not exceed the torque of 50 Ncm. Insert it in the implant seat and complete the placement of the implant, so that the laser treated portion is completely inserted in the bone tissue.



Insertion with micromotor



Remove the implant from the touch&go holder using the Beak driver.

Insert the implant, keeping below 15 rpm; do not exceed the torque of 50 N•cm.

Tightening of the cover screw



Remove the cover screw from the upper part of the touch&go holder, using a spanner with the Microesam tip. After having cleaned the implant seat, tighten the cover screw, with a maximum torque of 15 Ncm.

Surgical organizers

way organizer

instruments not included

31351

PPSU

SI

All instruments of way lines, including way Slim and way Short, are hosted in a unique tray, organized according to a logical and intuitive path.


Main characteristics:


- **functional:** thanks to the hinged lid, it is possible to incline it and to easily take the instruments;
- **rational:** in the tray there are the rotary instruments and the inserts; on bottom of the tray, the wrenches are accommodated;
- **safe:** once closed, the tray remains blocked to avoid any movement of the instruments;
- steam sterilisable up to 134°C

Supplied with x-ray template.



Arrangement of instruments inside way organizer







way Slim organizer

In a single organizer, all of the surgical and restorative instruments for way Slim implants. Supplied with x-ray template.

instruments not included

30665



way Short organizer


In a single organizer, all of the surgical and restorative instruments for way Short implants. Supplied with x-ray template.

instruments not included

30663

PPSU

SI



Drills

The visual references present on the drills allow you to evaluate the depth drilled based on the length of the implant chosen.

The drilling phases must be carried out with an up and down movement, without exceeding the maximum speed indicated in each phase of the protocol. The use of the Drill Controller for the twist drills and the Stop for the final drills facilitates the perforation.

Do not use drills which result as damaged, are not sharp or which have been used for more than 4 applications in order to reduce risks of overheating and bone trauma which may compromise the osseointegration process.

Drill extension

| | |
|------|-------|
| | 28450 |
| Inox | |

To be used with rotating instruments in order to easily reach the intervention regions between two dental elements.




Round bur

| | |
|------|-------|
| | 15917 |
| Inox | WC |

To be used as an alternative to the lance drill or to level any small unevenness on the bone crest.



Spade drill



| | |
|-------|-------|
| short | 30659 |
| long | 30661 |


Inox

This creates a niche on the cortical bone for the subsequent drills. It creates a precise entrance point thanks to its perfect centering

- Do not sink the instrument up to the final length of the implant to be inserted; use the reference notches to always maintain a margin of at least 2 mm between the depth of the drill and that of the implant site.



Indicator pin



| | | |
|-------|--|-------|
| ø 2.1 | | 15776 |
| ø 2.5 | | 13207 |

Ti

PU

Inserted into the implant site being created, it indicates axis and depth thanks to the notch, as shown in the side diagram.



Enlarger drill ø 3


| |
|-------|
| 25552 |
|-------|

Inox

To be used for way Slim implants in case of D1 bone; it should be sunk up to the mark, regardless of the implant length.



Twist drill




| | | | |
|---|-------|-------|-------|
| □ | ø 2.1 | short | 31355 |
| □ | ø 2.1 | long | 31254 |
| ■ | ø 2.5 | short | 31357 |
| ■ | ø 2.5 | long | 31256 |

Inox



This prepares the implant site based on the length of the chosen implant. The measurements indicated by the notches do not include the tip of the drill, about 0.7 mm. It is therefore advisable to consider this difference when planning the perforation phases. The drills are to be matched only with the dedicated stops, shown in the current catalogue; do not use other stops, as an implant site with wrong dimensions could thus be created, with serious risks for the patient.



Stop for twist drill



| | Implant length | | | | | | | | |
|-------|----------------|-------|-------|-------|-------|-------|-------|-------|-------|
| | 5 | 6.5 | 8 | 9 | 10 | 11 | 12 | 13 | 15 |
| short | 31259 | 31262 | 31265 | 31268 | 31271 | 31274 | 31277 | 31280 | 31283 |

TI



| | Implant length | | | | | | | | |
|------|----------------|-------|-------|-------|-------|-------|-------|-------|-------|
| | 5 | 6.5 | 8 | 9 | 10 | 11 | 12 | 13 | 15 |
| long | 31286 | 31289 | 31292 | 31295 | 31298 | 31301 | 31304 | 31307 | 31310 |

TI







Designed to reach the planned depth with the twist drills in a precise, safe and controlled manner, safeguarding the anatomical respect zones. Stops dedicated to long drills are distinguished by letter "L".

After having inserted the stop, always verify that the free part of the drill corresponds to the planned length.

After drilling, remove the stop, paying attention to the cutting edge of the drill.

! They must be exclusively used with the twist drills shown on the current catalogue, otherwise serious risks can be caused to the patient

Final drill


| | 5 | 6.5 | 8 | 9 | 10 | 11 | 12 | 13 | 15 |
|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|
|  | | | | 29643 | 31220 | 29646 | 29649 | 29652 | 29655 |
|  | 19740 | 19731 | 29658 | 29661 | 29664 | 29667 | 29670 | 29673 | 29676 |
|  | 19737 | 19728 | 29679 | 29682 | 29685 | 29688 | 29691 | 29694 | 29697 |
|  | 19734 | 19725 | | | | | | | |

 Inox

It allows you to complete the implant site with widening suitable for the dimensions of the implant.


The nominal depth of the drill (tip included) corresponds with the notch where the working part of the drill finishes.

The inferior limit of the coloured band corresponds to a depth increase of 1 mm; the upper limit corresponds to a depth increase of 3 mm.



The final drills L.5 and 6.5 mm come equipped with integrated drill stop, beyond which you must not descend.

On one side the stop has been milled with a longer cutting edge making it possible to level the bone crest.




Stop for final drill

| | |
|---|-------|
|  | 16549 |
|  | 10953 |
|  | 10954 |

 Inox

It makes it possible to control the perforation depth during the last phase of drilling, reducing risks to the anatomical respect zones.

The final drill stops must not be used with 5 and 6.5 mm drills as the stop is integrated into the drill itself.




Insert the stop on the drill from above, following the indicated direction. Fix the drill to the handpiece. Press the device in the part indicated by the arrows and push until it clicks into position.



Bone Profiler organizer

It hosts all the instruments Bone Profiler, besides the spanners and inserts necessary to the restoration with Mua abutments, thus avoiding to sterilise the complete surgical kit.

| | |
|------|-------|
| | 31541 |
| PPSU | SI |



Bone Profiler

| | |
|-----|-------|
| H 1 | 30720 |
| H 3 | 30722 |

Inox

It allows to level the bone crest, thus creating the necessary space to correctly place the Mua abutment.

It must be used in combination with the guide, into which it must be correctly inserted, before being started. Maximum speed of use: 200 rpm. Choose the instrument, based on the height of the Mua abutment to be used.





Bone Profiler guide

| | |
|----|-------|
| | 30727 |
| Ti | |


It keeps the bone profiler in the correct axis during its use, thus protecting the implant connection. Before using the bone profiler, screw the guide into the implant, tightening it with the Microesam instruments at 15 Ncm.



Protocol of Bone Profiler





After placement of the implant, tighten the guide for bone profiler into the implant seat with the Microesam instruments at 15 Ncm.



Choose the correct bone profiler, based on the height of the Mua abutment. Insert the bone profiler into the guide and, **only after positioning**, start it at the maximum speed of 200 rpm.

In case of an implant placement particularly under the crest, pay attention to the soft tissues; if necessary, open a flap to avoid damage.



After having leveled the bone crest, remove the guide and place the Mua abutment.

Spanners and inserts



To be used for handling the implants and prosthetic components.

All inserts can be used alone or in combination with the screwdriver, the ratchet wrench or the torque wrench; in the latter cases, verify that the matching between the two devices is correct.

The drivers are to be inserted on the handpiece to handle the various devices easily and quickly; ensure that they are effectively retained. A maximum speed of 15rpm is advisable.

For the tightening of the prosthetic components, always use a controlled torque wrench, as the use of the screwdriver or of the ratchet wrench can easily lead to excessive torque. When using spanners and inserts, it is important to avoid lateral bendings, which may cause the instrument break or the damage of the handled components.

Screwdriver



| | | |
|---|---|-------|
| | | 14242 |
| <i>o-ring</i> | <i>pack 3pcs</i> | 21143 |
|  |  | |

It allows you to use the various inserts manually, giving you the utmost perception and sensitivity in your handling.

You will feel a click when the insert connects with the screwdriver, indicating that insertion has taken place correctly.




Universal screwdriver



| | | |
|---|---|-------|
| | | 28641 |
| <i>o-ring</i> | <i>pack 3pcs</i> | 21143 |
|  |  | |

Matched with the inserts, it allows an easy handling during the implant insertion, thus guaranteeing an excellent control of direction.


Due to the high torque values it can easily reach, it must not be used for the tightening of the prosthetic components.



Insert extension

| | | |
|---|---|-------|
| | | 21126 |
| <i>o-ring</i> | <i>pack 3pcs</i> | 21144 |
|  |  | |

To be used with the inserts in order to easily reach the intervention region between two dental elements.



Newton screwdriver

| | | |
|-----------|-----------|-------|
| | | 26870 |
| lubricant | | 17002 |
| o-ring | pack 3pcs | 21143 |
| Inox | PTFE | PPS |

Coupled with the inserts, it is used to screw and unscrew the various devices easily and quickly, according to two modes: ratchet (no predefined torque limit) or dynamometric (calibrated torque).

Once the suitable insert has been selected, insert it into the head of the wrench and verify that the hexagonal profiles of the two devices are properly coupled. Rotate it until the hexagonal profiles of the two devices are coupled and push the insert downwards.

The arrow 'IN' on the wrench head shows the position of the wrench when tightening; by turning the device over, the arrow "OUT" is used when loosening screws. Pay attention to the position of the wrench during use, so that the rotation axis coincides with the axis of the handled device.

To select the torque limit, rotate the handle until the desired value; the adjustment key supplied in the package allows you to switch faster from one value to the other.




The selection of the torque must always be performed, during screwing of the handle clockwise; therefore, in order to adjust to an inferior torque to that set, it is important to unscrew the handle by two turns under the desired torque value, then screw again to the required value.

Tightening in torque mode

To tighten with a pre-set value of torque, rotate the handle until it is positioned exactly in line with the required value, then move the wrench in the direction indicated by the arrow; once the torque value has been reached, the wrench spins freely.

Tightening in ratchet mode


To use the wrench without a pre-set value of torque, rotate the handle until the writing "R".



Implant placement

► way Mix and way Extra

W-Start screwdriver




| | |
|---------------|-------|
| | 31507 |
| <i>o-ring</i> | 19034 |



To remove the implant from touch&go holder and insert it for some threads into the implant site manually. It differs from the W-Fix insert for the presence of the o-ring and for the fact that it cannot be used with the Newton torque wrench.




W-Start driver




| | |
|---------------|-------|
| short | 31510 |
| long | 31513 |
| <i>o-ring</i> | 19034 |



To remove the implant from touch&go holder and insert it for some threads into the implant site with the micromotor.



W-Fix insert



| | |
|-------|-------|
| short | 31516 |
| long | 31519 |




To complete the implant insertion into the implant site, together with the Newton torque wrench.



► way Slim

The instruments for way Slim are characterized by a black notch.

W-Start screwdriver ø 3




| | |
|---------------|-------|
| | 23491 |
| <i>o-ring</i> | 23485 |



To remove the implant from touch&go holder and insert it for some threads into the implant site manually. It differs from the W-fix insert for the presence of the o-ring and for the fact that it cannot be used with the Newton torque wrench.



W-Start driver ø 3




| | |
|---------------|-------|
| short | 23481 |
| long | 23484 |
| <i>o-ring</i> | 23485 |



To remove the implant from touch&go holder and insert it for some threads into the implant site with the micromotor.



W-Fix insert ø 3



| | |
|-------|-------|
| short | 23486 |
| long | 23487 |




To complete the implant insertion into the implant site, together with the Newton torque wrench.



► way Rock and way Short

Beak insert




| | |
|---------------|-------|
| short | 21207 |
| long | 21208 |
| <i>o-ring</i> | 23025 |


Inox

SI

To remove the implant from touch&go holder and insert it for some threads into the implant site



Beak driver



| | |
|---------------|-------|
| short | 21209 |
| long | 21097 |
| <i>o-ring</i> | 23025 |

Inox


SI

To remove the implant from touch&go holder and insert it into the implant site with the micromotor.



Tightening prosthetic components


Microesam insert




| | |
|-------|-------|
| short | 11655 |
| long | 11656 |



To be also used with the majority of the prosthetic components of way implants.




Microesam driver




| | |
|-------|-------|
| short | 11657 |
| long | 11658 |



To be also used with the majority of the prosthetic components of way implants.



Stepper insert



| | |
|-------|-------|
| short | 10473 |
| long | 10474 |



To handle the Mua straight abutment (way Mix) and the Reflect abutment (way Rock).



Equator insert

| | |
|----------------------|-------|
| | 26496 |
| holder (replacement) | 26497 |



To handle the Equator abutments.




Insertion-extractor tool Equator

| | |
|--|-------|
| | 26868 |
|--|-------|



To insert and remove the caps of the Equator system. Autoclavable.




Surgical planning

X-ray template

| | |
|----------------------|-------|
| way Mix and way Rock | 19054 |
| way Extra | 19053 |
| way Slim | 19465 |
| way Short | 20599 |


It shows all the implant sizes, according to the following scales:
1:1 Computerized Tomography (CT); 1,1:1 Endoral radiography
1,25:1 Orthopantomogram (OPG)



Depth probe

| | |
|----|-------|
| | 29993 |
| Ti | |

Ideal instrument to verify the depth of the osteotomy; the various sizes are also reported on the shank to facilitate the reading.



Equipment for implantology



Thanks to a detailed design, safety of use and high level materials, the W&H equipment satisfies the required standards, thus guaranteeing maximum precision in oral surgery.

The Implantmed motors guarantee high reliability; the contra angles are characterized by strength and effective ergonomics.


For those who use the ultrasound technology, Piezomed ensures maximum precision, without damaging the soft tissues. Thanks to Osstell Beacon, it is easy and fast to verify osseointegration and implant stability.



way|slim




way|miX



way|extra



way | miX




Ø 3.4

L

| | |
|----|-------|
| 9 | 29216 |
| 10 | 31217 |
| 11 | 29219 |
| 12 | 29222 |
| 13 | 29225 |
| 15 | 29228 |




1:1




Ø 3.8

L


| | |
|-----|-------|
| 6.5 | 31223 |
| 8 | 29231 |
| 9 | 29234 |
| 10 | 29237 |
| 11 | 29240 |
| 12 | 29243 |
| 13 | 29246 |
| 15 | 29249 |




1:1



| $\varnothing 4.5$ | |
|-------------------|-------|
| L | |
| 6.5 | 31226 |
| 8 | 29252 |
| 9 | 29255 |
| 10 | 29258 |
| 11 | 29261 |
| 12 | 29264 |
| 13 | 29267 |
| 15 | 29270 |



way|extra






ø 3.4

| L | |
|----|-------|
| 9 | 27281 |
| 11 | 27283 |
| 13 | 27285 |
| 15 | 27287 |



Ti






ø 3.8

| L | |
|----|-------|
| 9 | 23164 |
| 11 | 23055 |
| 13 | 23058 |
| 15 | 23061 |



Ti

































The implants way Extra d. 4.5 prior to 2020 (code 23167 - 23064 - 23067 - 23070) presented a different connection, equal to the implants way Milano d 4.5 and 5.5 mm.

! As the external morphology has not been modified, if both the new and previous version are present in the clinic, great attention must be paid to use the correct prosthetic components.

Prosthetic components way Mix and way Extra

Use way Milano restorative line for way Extra implants



| Healing | |  |  |  |  | | |
|-----------------------|--|---|---|--|--|---|---|
| | | Fit fixing screw | Finxing screw h 0.6 | Healing abutment | Wide healing abutment | | |
| | | single elements | | Bridges | | Structures | |
| Impression | |  |  |  |  |  |  |
| | | Pick-up coping | Fine Pick-up coping | Basic coping | Pick-up coping | Basic coping | Overdenture coping |
| Temporary restoration | |  |  |  | | | |
| | | Temporary abutment Single-Temp | Temporary abutment Esthetic-Temp | Temporary abutment Multi-Temp | | | |
| Cemented | |  |  |  |  |  | |
| | | Precision abutment | Custom abutment | Roller abutment | Elpy abutment | Moncone a finire | |
| Screwed | |  |  |  |  |  |  |
| | | Fusion abutment | Fusion CoCr abutment | Abutment for bar | Fusion CoCr abutment | Toronto abutment | Abutment for bar |
| | | | |  | |  |  |
| | | | | Mua abutment | | Mua abutment | Equator abutment |

Fixing screw is always supplied with the prosthetic components; this screw is to be used for the definitive fixing only.


► Management of soft tissues

Fit screw

Supplied with the implant.

It is ideal for cases with thin gingival biotype, as it closes flush with the implant seat, thus avoiding to hinder the healing of the soft tissues.

| 31335 | | | |
|-------|---|---------|------------|
| Ti |  | 15 N•cm | pack 3 pcs |



Cover screw h 0.6

Can be bought separately.

To be used in cases of adequate mucous thickness. The slight protrusion as compared with the implant seat prevents the growth of bone tissue over the screw itself, thus avoiding the subsequent difficulties of removal, in case the implant has been positioned below the crest.

| 31336 | | | |
|-------|---|---------|------------|
| Ti |  | 15 N•cm | pack 3 pcs |



Healing abutment

To guide the healing of periimplant soft tissues.

| H | T | T |
|---|-------|-------|
| 2 | 30191 | 29441 |
| 3 | 31342 | 31345 |
| 4 | 30193 | 29444 |
| 6 | 30195 | 29447 |
| 8 | 30197 | 29450 |

Ti  15 N•cm




Wide healing abutment

In cases of high mucous trait.

| Ø | H | T | T |
|-----|---|-------|---|
| 4.5 | 4 | 30199 | |
| 4.5 | 6 | 30201 | |
| 5.5 | 4 | 30203 | |
| 5.5 | 6 | 30205 | |
| 6.5 | 4 | 29459 | |
| 6.5 | 6 | 29462 | |

Ti  15 N•cm




Temporary abutment Single-Temp

For single temporary elements.

| H | T | T |
|-------------|-------|-------|
| 1 | 30207 | 29465 |
| 3 | 30209 | 29468 |
| 5 | 30211 | 29471 |
| short screw | | 15833 |
| long screw | | 17227 |

Ti 15 N·cm




Temporary abutment Multi-Temp

For temporary restoration on multiple elements.

| H | T | T |
|-------------|-------|-------|
| 1 | 30213 | 29474 |
| 3 | 30215 | 29477 |
| 5 | 30217 | 29480 |
| short screw | | 15833 |
| long screw | | 17227 |

Ti 15 N·cm




Temporary abutment Esthetic-Temp

Ideal for esthetic areas.

| H | T | T |
|-------------|-------|-------|
| 1 | 31227 | 29483 |
| 3 | 31228 | 29486 |
| short screw | | 15833 |
| long screw | | 17227 |

PEEK Ti 15 N·cm




► Impression

Pick-up coping

For custom impression tray.

| H | T | T |
|------------------|-------|-------|
| 3 | 30219 | 29489 |
| screw | | 31544 |
| 7 | 30221 | 29492 |
| screw | | 30869 |
| screw extra long | | 29352 |


Ti 15 N·cm



Fine Pick-up coping

For custom impression tray, in cases of very near adjacent teeth.


| | T | T |
|------------------|-------|---------|
| | 31339 | |
| screw | 30869 | |
| screw extra long | 29352 | |
| | Ti | 15 N•cm |



Overdenture coping

For custom impression tray, it facilitates the impression taking in case of structures thanks to the absence of antirotational index.


| | T | T |
|-------|-------|---------|
| | 30223 | 29564 |
| screw | 30869 | |
| | Ti | 15 N•cm |



Basic coping

For standard impression tray.


| | T | T |
|-------|-------|---------|
| | 31532 | 31535 |
| screw | 17225 | |
| | Ti | 15 N•cm |



Cap for Basic coping

To be used with Basic coping to increase the accuracy of the impression, in specific cases of disparallelism.

| | |
|-----|-------------|
| | 16390 |
| POM | pack 10 pcs |



Analog

It reproduces the position of the implant in the plaster model.

| | T | T |
|--|-------|---|
| | 31247 | |
| | Ti | |




► Definitive restoration

Precision abutment

The versatile conformation makes them suitable for a wide variety of restorative solutions.


| H | T | T |
|-------|---------------|---------|
| | straight | |
| 1 | 30225 | 29498 |
| 3 | 30227 | 29501 |
| 5 | 30229 | 29504 |
| | angulated 15° | |
| 1 | 30231 | 29507 |
| 3 | 30233 | 29510 |
| 5 | 30235 | 29513 |
| | angulated 25° | |
| | 30237 | 29516 |
| | 30239 | 29519 |
| | 30241 | 29522 |
| screw | 15833 | |
| Ti | | 25 N·cm |



Custom abutment

The possibility to shape a personalized abutment allows a great flexibility of use, while at the same time maintaining the precision of the industrial coupling.

| T | T |
|-------|---------|
| | 30243 |
| screw | 15833 |
| Ti | 25 N·cm |



Milled Custom abutment

The specific morphology optimizes the milling operations; in some techniques, it can also be used as coping.

| T | T |
|-------|---------|
| | 31337 |
| screw | 15833 |
| Ti | 25 N·cm |



Elpy abutment


It allows to start the preparation of the crown as near as possible to the implant platform.

It cannot be used with way Extra implants.

The Small version is dedicated to implants d. 3.4 and 3.8 mm.

The Large version is to be used exclusively on implants d. 4.5 mm


| | T | T |
|---------------|---------|-------|
| straight | 29568 | 29552 |
| angulated 15° | 29453 | 29555 |
| angulated 25° | 29456 | 29558 |
| screw | 15833 | |
| Ti | 25 N•cm | |



Roller abutment

The particular design of the collar enhances and respects the anatomical profile of the soft tissues on the esthetic zones, while the pre-shaped margins facilitate laboratory works, reducing production times. The profile of the collar is oriented according to the position of the implant seat, which can be planned by using the W-fix insert.

| H | T | T |
|-------|---------------|-------|
| | straight | |
| 1 | 30245 | 29525 |
| 3 | 30247 | 29528 |
| 5 | 30249 | 29531 |
| | angulated 15° | |
| 1 | 30251 | 29534 |
| 3 | 30253 | 29537 |
| 5 | 30255 | 29540 |
| | angulated 25° | |
| 1 | 30257 | 29543 |
| 3 | 30259 | 29546 |
| 5 | 30261 | 29549 |
| screw | 15833 | |
| Ti | 25 N•cm | |



Shoulderless abutment

Specific for the vertical preparation technique


| | T | T |
|-------|---------|-------|
| | 30263 | 29561 |
| screw | 15833 | |
| Ti | 25 N•cm | |



Fusion abutment CoCr

For the creation of restorations with the casting technique.


| H | T | T |
|-------|--------------|---------|
| | non rotating | |
| 1 | 31319 | |
| 3 | 31322 | |
| | rotating | |
| 1 | 31313 | |
| 3 | 31316 | |
| screw | 15833 | |
| PMMA | Ti | CoCr |
| | | 25 N·cm |



Fusion abutment

For the creation of restorations with the casting technique.


| H | T | T |
|-------|-------|---------|
| 1 | 31229 | |
| 3 | 31230 | |
| screw | 15833 | |
| PMMA | Ti | Au lega |
| | | 25 N·cm |



Abutment for bar

Specific for overdenture restoration with bars.


| H | T | T |
|----------|-------|---|
| 1 | 31231 | |
| 3 | 31232 | |
| screw | 15835 | |
| castable | PMMA | 4 N·cm in lab 25 N·cm definitive element |
| base | Ti | W-FIX 25 N·cm |



Toronto abutment

Ideal for the creation of the toronto-bridge, ensuring perfect structure fit in the oral cavity.

| H | T | T |
|-----------------|-------|---|
| 1 | 31233 | |
| 3 | 31234 | |
| screw | 15835 | |
| castable spacer | 20210 | |
| cylinder | Ti | PMMA PEEK 4 N·cm in lab 25 N·cm definitive element |
| base | Ti | W-FIX 25 N·cm |



Equator abutment

To be used with the dedicated caps, to anchor the removable prostheses.

| H | T | T |
|---|-------|---|
| 1 | 31235 | |
| 2 | 31236 | |
| 3 | 31237 | |
| 4 | 31238 | |
| 5 | 31239 | |
| 6 | 31240 | |

Ti 25 N·cm



Equator cap

It allows to correct disparallelism up to 25°.

Each pack contains:

1 container for caps in titanium

1 black cap for lab use

1 protective disk

4 retentive caps (1 for each retention grade)

| |
|-------|
| 26861 |
|-------|

Ti PA EVA



Smartbox kit

It allows to correct disparallelism up to 50°.

Each pack contains:

1 container with cap for lab

1 pink protective disk

4 retentive caps (1 for each retention grade)

| |
|-------|
| 27723 |
|-------|

Ti PA EVA



Spare containers

| | | |
|------------|-------|---------------------|
| | | |
| Inox | Ti | Smartbox (with cap) |
| 24088 | 24089 | 27724 |
| pack 2 pcs | | pack 1 pc |

spare caps

| | | | | | |
|----------------|-----------------|------------|--------------|--------------|----------------|
| | | | | | |
| white/clear | yellow | pink | violet | black | black Smartbox |
| Standard 1800g | Extra-soft 600g | Soft 1200g | Strong 2700g | Only for lab | Only for lab |
| 26864 | 26863 | 26865 | 26862 | 24087 | 27725 |
| pack 4 pcs | | | | | |

► Mua

Created for the total fixed rehabilitation on distally tilted implants, the MUA components allow the emerging parts of inclined implants in posterior sectors to be parallel, thus simplifying the prosthesis positioning, fitting and fixing. Available in the straight version and angulated at 17° and 27°, the range allows to choose between two collar heights (1 and 3mm).

The use of MUA abutments requires dedicated prosthetic components. In order to correctly place the MUA abutments, in many cases it is necessary to use the bone profiler to level the bone crest and create the necessary space (see page 28).



Mua straight abutment

For easy handling, there is an accessory in peek in the pack. Once the straight abutment is positioned, simply bend and remove the accessory, pulling it out; then, fix the abutment with the Stepper insert at the indicated torque.


| H | T | T |
|----|-------|---------|
| 1 | 31241 | |
| 3 | 31242 | |
| Ti | PEEK | 25 N·cm |



Mua angulated abutment

It includes a titanium pre-mounted accessory, which facilitates the positioning and allows to verify the direction of the prosthetic axis. Once tightened the angulated abutment with the Microesam insert at the indicated torque, remove the accessory by unscrewing it for a few rounds.

| H | T | T |
|-------|---------------|---------|
| | angulated 17° | |
| 1 | 31243 | |
| 3 | 31244 | |
| | angulated 27° | |
| 1 | 31245 | |
| 3 | 31246 | |
| screw | 25868 | |
| Ti | | 25 N·cm |



Mua healing abutment

Used during the healing phase of soft tissues.

| T | T |
|-------|---------|
| 25848 | |
| Ti | 15 N·cm |



Mua analog

It recreates the position of the implant, on which the Mua abutment has been fixed.


| | | |
|----|-------|---|
| | T | T |
| | 25851 | |
| Ti | | |



Mua Pick-up coping

For custom impression tray.


| | | |
|-------|---------|---|
| | T | T |
| | 29970 | |
| screw | 25974 | |
| Ti | | |
| | 15 N·cm | |



Mua App abutment

It may act as temporary abutment, coping or as definitive solution.


| | | |
|-------------|---------|---|
| | T | T |
| | 25854 | |
| short screw | 25865 | |
| long screw | 25974 | |
| Ti | | |
| | 15 N·cm | |



Mua App accessories

It consists of castable part and spacer, to create a definitive restoration with the App abutment.


| | | |
|------|-------|---|
| | T | T |
| | 26871 | |
| PMMA | PEEK | |




Mua castable

To create the definitive restoration.


| | | |
|-------|----------------------------|---|
| | T | T |
| | 25862 | |
| screw | 25865 | |
| PMMA | Ti | |
| | 4 N·cm in lab | |
| | 15 N·cm definitive element | |









way | slim



| | ø 3 | Ti |
|----|-------|----|
| L | | |
| 9 | 30479 | |
| 11 | 25405 | |
| 13 | 25406 | |
| 15 | 25407 | |



way Slim prosthetic components

| | | | |
|------------------------|--|---|--|
| Healing | Healing abutment  | | |
| | Single elements | Bridges | Structures |
| Impression | Pick-up coping  | | |
| Temporary restoration | Temporary abutment Single-Temp  | Temporary abutment Multi-Temp  | |
| Definitive restoration | Precision abutment  | | Equator abutment  |

Fixing screw is always supplied with the prosthetic components; this screw is to be used for the definitive fixing only.

► Management of soft tissues

Cover screw


Supplied with the implant.



Healing abutment

To guide the healing of periimplant soft tissues.


| H | |
|---|-------|
| 2 | 23475 |
| 4 | 23476 |
| 6 | 23477 |



Temporary abutment Single-Temp

For single temporary elements.


| H | |
|-------------|-------|
| 1 | 23452 |
| 3 | 23455 |
| screw short | 15833 |
| screw long | 17227 |



Temporary abutment Multi-Temp

For temporary restoration on multiple elements.

| H | |
|-------------|-------|
| 1 | 20561 |
| 3 | 20564 |
| screw short | 15833 |
| screw long | 17227 |




► Impression

Pick-up coping

For custom impression tray.

| | |
|-------|-------|
| | 23478 |
| screw | 17227 |



Analog

It reproduces the position of the implant on the plaster model.

| | |
|--|-------|
| | 25408 |
|--|-------|




► Definitive restoration

Precision abutment

Suitable for a wide variety of restorative solutions.

| H | |
|---------------|-------|
| straight | |
| 1 | 23458 |
| 2 | 23461 |
| 4 | 23464 |
| angulated 15° | |
| 1 | 23467 |
| 2 | 23470 |
| 4 | 23473 |
| screw | 15833 |

Ti   25 N•cm



Equator abutment

To be used for rehabilitation with removable prostheses, together with the caps and the dedicated accessories (see page 49).

| H | |
|---|-------|
| 1 | 26472 |
| 2 | 26475 |
| 3 | 26478 |
| 5 | 26481 |

Ti   25 N•cm








way|rock





way|short




| Ø 3.4 | |
|--------------|-------|
| platform 4.8 | |
| L | |
| 11 | 30585 |
| 12 | 30588 |
| 13 | 30591 |
| 15 | 30594 |



| Ø 3.8 | |
|--------------|-------|
| platform 4.8 | |
| L | |
| 8 | 30597 |
| 9 | 30600 |
| 10 | 30603 |
| 11 | 30606 |
| 12 | 30609 |
| 13 | 30612 |
| 15 | 30615 |








ø 4.5
platform 4.8





| L | |
|----|-------|
| 8 | 30618 |
| 9 | 30621 |
| 10 | 30624 |
| 11 | 30627 |
| 12 | 30630 |
| 13 | 30633 |
| 15 | 30636 |




ø 4.5
platform 6.5




| L | |
|----|-------|
| 8 | 30639 |
| 9 | 30642 |
| 10 | 30645 |
| 11 | 30648 |
| 12 | 30651 |
| 13 | 30654 |
| 15 | 30657 |





way | short



| |
|-------------------|
| \varnothing 3.8 |
| platform 4.8 |
| L |
| 5 27274 |
| 6.5 27277 |



| |
|-------------------|
| \varnothing 4.5 |
| platform 4.8 |
| L |
| 5 19788 |
| 6.5 19782 |




















| |
|-------------------|
| \varnothing 5.5 |
| platform 6.5 |
| L |
| 5 19785 |
| 6.5 19779 |



► way Rock and way Short prosthetic components

Use way Rock restorative line for way Short implants

| | | | | |
|------------------------|------------|--|---------|------------|
| | Healing |  Healing abutment | | |
| | | Single elements | Bridges | Structures |
| | Impression |  Pick-up coping  Basic coping | | |
| | |  Temporary abutment Single-Temp  Temporary abutment Multi-Temp | | |
| Definitive restoration | Cemented |  Precision abutment  Rotative abutment  Custom abutment  Slender abutment | | |
| | Screwed |  Fusion abutment  Fusion CoCr abutment  Abutment for bar  Fusion CoCr abutment  Toronto abutment  Abutment for bar  Reflect abutment  Equator abutment | | |

Fixing screw is always supplied with the prosthetic components; this screw is to be used for the definitive fixing only.

► Management of soft tissues

Cover screw

Supplied with the implant.



Healing abutment

To guide the healing of periimplant soft tissues.

| H | T | T |
|---|-------|-------|
| 2 | 31358 | 31361 |
| 3 | 31359 | 31362 |
| 5 | 31360 | 31363 |




Single-Temp abutment

For single temporary elements.

| H | T | T |
|-------------|-------|-------|
| in PEEK | | |
| | 31364 | 31365 |
| short screw | 14026 | |
| long screw | 18735 | |



| H | T | T |
|-------------|-------|-------|
| in TITANIUM | | |
| | 31366 | 31367 |
| short screw | 14026 | |
| long screw | 18735 | |




Multi-Temp abutment

For temporary restoration on multiple elements.

| H | T | T |
|-------------|-------|-------|
| in PEEK | | |
| | 31368 | 31369 |
| short screw | 14026 | |
| long screw | 18735 | |



| H | T | T |
|-------------|-------|-------|
| in TITANIUM | | |
| | 31370 | 31371 |
| short screw | 14026 | |
| long screw | 18735 | |




► Impression

Pick-up coping

For custom impression tray.

| H | T | T |
|-------|-------|-------|
| 3 | 31372 | 31374 |
| screw | 20066 | |
| 6 | 31373 | 31375 |
| screw | 18735 | |

Ti   15 N·cm




Basic coping

For standard impression tray.

| H | T | T |
|-------|-------|-------|
| | 31376 | 31377 |
| screw | 20066 | |

Ti   15 N·cm




Cap for Basic coping

To be used with Basic coping to increase the accuracy of the impression, in specific cases of disparallelism.

| | |
|--|-------|
| | 16390 |
|--|-------|

POM *pack 10 pcs*



Analog

It reproduces the position of the implant in the plaster model.

| H | T | T |
|---|-------|-------|
| | 31378 | 31379 |

Ti




► Definitive restoration

Precision straight abutment

Versatile conformation; suitable for a wide variety of restorative solutions.

| H | T |
|-------|-------|
| 1 | 31380 |
| 2 | 31381 |
| 3 | 31382 |
| screw | 14026 |

Ti   35 N·cm




Rotative abutment

Made from two distinct components, they allow for a rotation of 360° of the coronal part of the abutment. Once the best position has been selected, the two components are fixed with the passing screw.

| H | T | T |
|-------|---------------|-------|
| | angulated 15° | |
| 1 | 31383 | 31387 |
| 3 | 31384 | 31388 |
| screw | | 18736 |
| | angulated 25° | |
| 1 | 31385 | |
| 3 | 31386 | |
| screw | | 18736 |

Ti   35 N·cm



Custom abutment

It allows for maximum personalization, while at the same time maintaining the precision of the industrial coupling.

| | T | T |
|-------|-------|-------|
| | 31389 | 31390 |
| screw | | 14026 |

Ti   35 N·cm



Slender abutment

It allows you to manage cases where the gingival tissue tends to have a very thin thickness.






| H | T | T |
|-------|---------------|-------|
| | straight | |
| | 31391 | 31394 |
| screw | | 14026 |
| | angulated 15° | |
| | 31392 | 31395 |
| | angulated 25° | |
| | 31393 | 31396 |
| screw | | 18740 |


Ti   35 N·cm



Fusion abutment Co Cr

For the creation of restorations with the casting technique.


| | T | T |
|--------------|---|---|
| non rotating |  |  |
| | 31328 | 31334 |
| rotating |  |  |
| | 31325 | 31331 |
| screw | 14026 | |
| PMMA | Ti | CoCr |
| |  | 35 N•cm |



Fusion abutment



For the creation of restorations with the casting technique.


| H | T | |
|-------|---|---------|
| 1 | 31397 | |
| screw | 14026 | |
| PMMA | Ti | Au lega |
| |  | 35 N•cm |



Abutment for bar

Specific for overdenture restoration with bars.


| H | T | T |
|----------|--|---|
| 1 | 31398 | 31400 |
| screw | 18738 | |
| 3 | 31399 | 31401 |
| screw | 18739 | |
| castable | PMMA |  |
| | 4 N•cm in lab 35 N•cm prosthetic artifact | |
| base | Ti |  |



Toronto abutment


Ideal for the creation of the toronto-bridge, ensuring perfect structure fit in the oral cavity.

| | T | |
|-----------------|---|---------|
| | 31402 | |
| screw | 14026 | |
| castable spacer | 20213 | |
| Ti | PMMA | PEEK |
| |  | 35 N•cm |



Compact abutment

Characterized by its mechanical solidity which comes from its full structure and the conical coupling



| | T | T |
|-----|---|---------|
| 4 | 31403 | 31406 |
| 5,5 | 31404 | 31407 |
| 7 | 31405 | |
| Ti |  | 35 N•cm |




Equator abutment

To be used for rehabilitation with removable prostheses, together with the caps and the dedicated accessories (see page 49)

| | T |
|---|-------|
| 1 | 31408 |
| 2 | 31409 |
| 3 | 31410 |
| 4 | 31411 |
| 5 | 31412 |
| 6 | 31413 |

Ti


35 N·cm




▶ Reflect


Reflect abutment

Designed for the creation of screw-retained restorations. Once tightened on the implant, it must not be removed; it requires dedicated components.

| | T | T |
|--|-------|-------|
| | 31414 | 31415 |

Ti


35 N·cm




Reflect Pick-up coping

For custom impression tray.

| | T | T |
|-------|-------|-------|
| | 31416 | 31417 |
| screw | | 20711 |

Ti


15 N·cm



Reflect analog

It reproduces the position of the implant with the Reflect abutment on the plaster.

| | T | T |
|--|-------|-------|
| | 31418 | 31419 |



Ti




Reflect Melty abutment

To create a definitive restoration on the Reflect abutment.

| | T | T |
|-------|-------|-------|
| | 31420 | 31421 |
| screw | | 18738 |

Ti
PMMA


4 N·cm in lab
35 N·cm prosthetic artifact



digitals

Digital evolution

The digital solutions proposed by Geass constitute an open, flexible and economical system that allows you to:

- acquire a complete digital flow;
- utilise innovative technologies that are adaptable to any IT system;
- implement your own digital structure with versatile equipment and instruments;
- define an effective surgery-laboratory work flow.

Geass is able to offer complete, professional service for your informed investment in the digital world, with specific consultancy, operational support in your surgery and in the lab, events and training courses and after-sales assistance on-site and remotely.

Digital impression taking

Intraoral scanner


Ease of use, fluid process of image acquisition, precision results: all this for an excellent workflow, accurate and fast.



Computer assisted surgery

Gedrive

Gedrive is the computer assisted surgery which allows you to carry out accurate and complete diagnoses, plan the optimal position of the implant and design the prosthetic rehabilitation while operating in total safety.



cad-cam personalized prosthesis

Performa

To overcome the limits of the current CAD-CAM productions on implants, Geass uses an innovative technology which integrates the advantages of milling with those of machine turning. Various solutions in several types of Zirconia, PMMA and laser melting complete the Performa offer.



Gedrive

Gedrive way organizer

instruments not included


28475



PPSU

SI

Gedrive Bluebox is an organizer designed to contain in a rational and functional manner:

- the Gedrive Start drills, to carry out guided surgery in the initial phases;
- the drills, wrenches and inserts of the Final Gedrive, to be assisted right up to implant placement.





Drill for pin

| | |
|------|-------|
| | 30437 |
| Inox | |

Allows you to create a seat for the fixing pins.; maximum speed: 500 rpm



Fixing pin

| | |
|----|-------|
| | 20341 |
| Ti | |

Allows you to fix the surgical guide.



Mucotome

| | |
|-------|-------|
| ø 3.8 | 23242 |
| ø 4.5 | 23245 |
| Inox | |

To incise and remove the soft tissues; maximum speed: 40 rpm.



Centering drill

| | | |
|-------|--|-------|
| ø 2.1 | | 23248 |
| ø 2.8 | | 23272 |
| Inox | | |

It creates the first osteotomy to facilitate the precise centering and positioning for the subsequent drills, thus levelling the bone crest at the same time, if necessary. Maximum speed: 400 rpm.






Twist drill


| | | | | | | | | |
|-------|--|-------|-------|-------|-------|-------|-------|-------|
| | | 8 | 9 | 10 | 11 | 12 | 13 | 15 |
| ø 2.1 | | 23251 | 23254 | 23257 | 23260 | 23263 | 23266 | 23269 |
| ø 2.8 | | 23275 | 23278 | 23281 | 23284 | 23287 | 23290 | 23293 |
| Inox | | | | | | | | |

For the initial preparation of implant site; the integrated stop guarantees more safety. Maximum speed: 400 rpm




Final drill




| | 8 | 9 | 10 | 11 | 12 | 13 | 15 |
|---|-------|-------|-------|-------|-------|-------|-------|
|  | | 23296 | | 23299 | 23302 | 23305 | 23308 |
|  | 23311 | 23314 | 23317 | 23320 | 23323 | 23326 | 23329 |
|  | 23430 | 23332 | 23335 | 23338 | 23341 | 23344 | 23347 |




The final drill allows you to complete the implant site with widening adequate to the dimensions of the implant; maximum speed: 300 rpm.



Enlarger drill




| | |
|---|-------|
|  | 23549 |
|  | 23552 |
|  | 23555 |






To be used in cases of D1 bone; maximum speed: 300 rpm.



W-Start screwdriver




| | | |
|---|---|---|
|  |  |  |
| 23368 | 31520 | |
| <i>o-ring. (3 pcs)</i> | 15928 | 15928 |




  

To remove the implant from touch&go holder and insert it for some threads into the implant site. It differs from the W-Fix insert for the presence of the o-ring and for the fact that it cannot be used with the Newton screwdriver.



W-Start driver




| | | |
|---|---|---|
|  |  |  |
| 23374 | 31526 | |
| <i>o-ring. (3 pcs)</i> | 15928 | 15929 |



  

To remove the implant from touch&go holder and insert it for some threads into the implant site. It differs from the W-Fix insert for the presence of the o-ring.




W-Fix insert

| | | |
|---|---|---|
|  |  |  |
| 23356 | 31523 | |

To be used with the screwdriver and the Newton torque wrench to complete the implant insertion into the implant site.



The other instruments hosted in the Geadrive organizer are in common with the traditional surgery (pages 24 to 35)





Contra angle adapter

| | |
|---|-------|
| | 29358 |
| <i>o-ring. (3 pcs)</i> | 21144 |
|   | |

To handle the mounters with the micromotor.



Mounter

| | |
|---|-------|
|    | |
| 31550 | 31531 |
|  | |

It allows to remove the implant from the touch&go holder and to place it into the implant site; it can be used with the screwdriver and the Newton screwdriver.

Do not exceed the torque of 50 Ncm. Left into the implant site, it helps to maintain the surgical guide in the correct position until the end of the intervention.




Mounter extractor

| | |
|---|-------|
| | 28305 |
|  | |

Screwed into the mounter instead of the screw, it allows to remove it, in case it remains blocked in the implant seat.



Sleeve for guide

| | |
|---|-------|
| ø 4.2 | 26511 |
| ø 5.2 | 26514 |
|  | |

Fixed on the surgical guide, it allows to guide the drills so that the osteotomy corresponds to the virtual planning of the treatment. The sleeve d. 4.2 is to be used with the implants d. 3.4 and 3.8 mm; the sleeve d. 5.2 is to be used with the implants d. 4.5 mm.



Sleeve for pin

| | |
|---|-------|
| | 28047 |
|  | |

Fixed on the surgical guide, it allows to guide the drill for pin.



Performa

Scanbody

To transfer the position of the implant from reality to the CAD software in three dimensions. They always have to be matched up with Geass library; the use of the matting spray is not required during scanning;

To tighten with Performa Torque at 4 N•cm.

As it is sterilizable, it can also be used for intraoral scanning; in this case use the Performa insert for fixing.




| way Mix - way Extra | way Slim | way Rock - way Short | | Mua | Reflect | |
|---------------------|-------------|----------------------|---------|-------|---------|---------|
| | | T p.4.8 | T p.6.5 | | T p.4.8 | T p.6.5 |
| short 31153 | short 25331 | 31435 | 31436 | 26551 | 32218 | 32219 |
| long 31154 | long 27251 | | | | | |
| | | | | | | |

Linker

Bases in titanium, on which it is possible to create CAD-CAM ceramic elements, especially suitable for cases of high aesthetic value. Characterised by knurling which facilitates retention of the cement. Colouring is yellow to reduce the metal reflection in transparency and to therefore improve the aesthetic outcome.

The height of the Linker is easily adaptable to the clinical situation, thanks to the pre-cut groove which, finding correspondence in the libraries, also facilitates the technician in the design of the prosthesis.




| | way Mix - way Extra | way Slim | way Rock - way Short | | Mua |
|------|---------------------|----------|----------------------|---------|------------|
| | | | T p.4.8 | T p.6.5 | |
| | H0 31156 | 29375 | 31441 | 31442 | H0 29088 |
| | H3 31157 | | | | |
| | H0 31158 | 29378 | 31443 | 31444 | H1,5 29089 |
| | H3 31159 | | | | |
| N•cm | 25 | 25 | 35 | 35 | 25 |




Linkers are available in two versions:

- V3: for way Mix, way Extra, way Rock and way Short; they have total height 10 mm and three pre-cutting grooves
- V2: for way Slim and Mua; they have total height 6 mm and a pre-cutting groove at 4 mm.

! Always verify the correct matching between the Linker to be used and the respective library



Linker V3



Linker V2

Base for wax-up

To be used in the laboratory to model the wax-up, then to be sent subsequently to Geass for digitalization.

| way Mix - way Extra | way Slim | way Rock - way Short | | Mua |
|---------------------|----------|---|---|-------|
| | |  p.4.8 |  p.6.5 | |
| 31160 | 25421 | 31439 | 31440 | 27308 |



pack 10 pcs



Ti-Base Sirona

| | way Mix - way Extra | way Slim | way Rock - way Short | |
|------|--|--|---|---|
| | | |  p.4.8 |  p.6.5 |
| | 31430 | 25632 | 31445 | 31446 |
| N•cm | 25  | 25  | 35  | 35  |



Digital analog

Specific for processing that derive from taking digital impressions, it ensures correct repositioning on the 3D printed model thanks to the presence of hexagonal sides, which also facilitate insertion.


The screw, included in the pack, ensures analog stability in the model and in many cases avoids the use of adhesive substances. The screw must be tightened with the Performa Torque tool.

Geass digital analogs are complete with the implant libraries necessary for virtual modelling through the main CAD softwares and for the creation of models using 3D printing.

| way Mix - way Extra | way Slim | way Rock - way Short | | Mua |
|---------------------|----------|---|---|-------|
| | |  b.4,8 |  b.6,5 | |
| 31155 | 28249 | 31437 | 31438 | 28261 |




Performa screwdriver




23918

Inox

To be used in the oral cavity to handle the scanbody on the implants.




Performa driver




28472

Inox

To be used in the oral cavity to handle Geass scanbody on the implants with micromotor.




Inclined hole insert




| | |
|-------|-------|
| short | 25449 |
| long | 25112 |

Inox

For tightening screws on tilted holes, used in CAD-CAM prosthesis. The definitive tightening of the screws on tilted holes is foreseen at 25 Ncm, except the screws on Mua which are to be tightened at 15 Ncm.



Inclined hole driver



| | |
|-------|-------|
| short | 25455 |
| long | 25452 |

Inox

For tightening screws on tilted holes, used in CAD-CAM prosthesis, with micromotor. The definitive tightening of the screws on tilted holes is foreseen at 25 Ncm, except the screws on Mua which are to be tightened at 15 Ncm.



Performa Torque

23788

Inox

POM

To be used exclusively in the lab to tighten the scanbody and the PMMA castables on the analogs, at a pre-defined torque of 4 Ncm.



advanced

Regenerative solutions

Easy Bone Management

EBM organizer

instruments not included 30377



For storing and organizing EBM drills.



Spherical EBM drill ø 5

29317



For sinus elevation with EBM technique, in compact bone.



Spherical EBM drill ø 7

29320



For sinus elevation with EBM technique, in soft bone.




Oval EBM drill ø 6

29323



To model the crestal portion of maxillary bones.



Sinus lift

Osteotome organizer

| | | |
|---------------------------------|----|-------|
| <i>instruments not included</i> | | 16499 |
| PPSU | PP | |

To store and organise the osteotomes simply.



Osteotome insert

| L | | | |
|----|-------|-------|-------|
| 9 | 11460 | 11466 | 11472 |
| 10 | 11461 | 11467 | 11473 |
| 11 | 11462 | 11468 | 11474 |
| 12 | 11463 | 11469 | 11475 |
| 13 | 11464 | 11470 | 11476 |
| 15 | 11465 | 11471 | 11477 |
| | Ti | | |

Specific for the minor sinus lift technique with way implants, it allows you to fracture the sinus cortical bone, raising it along with the membrane.



Osteotome handle

| | |
|------|-------|
| | 11459 |
| Inox | |

The insert suitable to the type of selected implant is to be inserted onto the hand piece.



Angled adapter for osteotome

| | |
|------|-------|
| | 16496 |
| Inox | |

Inserted onto the hand piece, it allows to reach the less easy areas.



Surgical mallet

| | |
|------|------------|
| | 2995.Y0.05 |
| Inox | |

It allows you to calibrate the force and concentrare it in a short instant, thus obtaining a precise fracture.



Endosteal elevator

| | |
|------|-------|
| n.1 | 11479 |
| n.2 | 11480 |
| n.3 | 11481 |
| Inox | |

Ideal to easily access the antral window, allowing you to carry out elevation of the sinus endosteum.




Split crest

Chisel double edge

| width | |
|-------|-------|
| 3 | 11485 |
| 6 | 11486 |
| 9 | 11487 |

Inox

Used to increase the transverse volume in greenstick osteotomy interventions.



Widener organizer

| | |
|---------------------------------|-------|
| <i>instruments not included</i> | 30545 |
|---------------------------------|-------|

PPSU

SI

Tray to effectively store the wideners and the handling instruments. The serigraphy allows to immediately identify the widener and the instruments to be used.




Widener

| No. | ø pilot hole | ø | ø | ø | ø | ø | ø |
|---------|--------------|------|------|-------|-------|-------|----------|
| | ø apical | L. 8 | L. 9 | L. 10 | L. 11 | L. 12 | L. 13-15 |
| 1 15702 | 2.0 | 2.20 | 2.25 | 2.30 | 2.35 | 2.45 | 2.50 |
| 2 15703 | 2.5 | 2.60 | 2.65 | 2.70 | 2.75 | 2.85 | 2.90 |
| 3 15704 | 2.8 | 2.95 | 3.00 | 3.10 | 3.15 | 3.25 | 3.30 |
| 4 15705 | 3.0 | 3.15 | 3.20 | 3.30 | 3.35 | 3.45 | 3.50 |
| 5 15706 | 3.5 | 3.60 | 3.70 | 3.75 | 3.85 | 3.90 | 4.00 |
| 6 15707 | 4.0 | 4.15 | 4.20 | 4.30 | 4.35 | 4.45 | 4.50 |

Ti

They allow you to gradually enlarge the crest, expanding the available bone and reducing surgical trauma. They increase the transverse volume in presence of thin edentulous ridges with suitable height.




Insert extension

| | |
|---------------|------------------------|
| | 21126 |
| <i>o-ring</i> | <i>pack 3pcs</i> 21144 |

Inox

SI

It allows you to use the wideners in less accessible areas or between two dental elements.



Bone regeneration

Tack holder

| | |
|--------------------|-------|
| tacks not included | 11342 |
| Ti | |

To conserve and organise the tacks tidily.



Tack inserter

| | |
|----|-------|
| | 11392 |
| Ti | Inox |


For the removal of the tack from the support and for positioning it *in situ* in conjunction with the mallet.



Tack

| | | |
|----|---|-----------|
| | L | |
| | 3 | 11369 |
| | 5 | 11370 |
| Ti | | pack 3pcs |

Used for fixing membranes.



SQ12 SQ17 organizer

| | |
|--------------------------|-------|
| instruments not included | 11341 |
| PPSU | Al |

It contains instruments and microscrews for fixing grafts, meshes and plates.



Drill SQ12 ø 1.0


| | |
|-------|-------|
| short | 11355 |
| long | 11356 |
| Inox | |

For creating an entrance on the cortical for the microscrews SQ12.

Manual spanner

| | |
|------|-------|
| SQ12 | 11344 |
| SQ17 | 11345 |

Inox



It provides excellent sensitivity when inserting the microscrew while applying adequate torque to it.

Contra-angle

| | |
|------|-------|
| SQ12 | 11346 |
| SQ17 | 11347 |

Inox




Extremely useful in areas which are difficult to reach thanks to the fitting on the hand piece.

Screwdriver

| | |
|------|-------|
| SQ12 | 11348 |
| SQ17 | 11349 |

Inox




To handle microscrews in the frontal region, allowing excellent directional control and of the torque applied.

Microscrew SQ12 ø 1.2

| | |
|----|-------|
| L | |
| 3 | 11371 |
| 4 | 11372 |
| 5 | 11373 |
| 6 | 11374 |
| 7 | 11375 |
| 8 | 11376 |
| 11 | 11377 |

Ti

pack 3pcs




Emergency microscrew
SQ12 ø 1.5

| | |
|---|-------|
| L | |
| 4 | 11378 |
| 6 | 11379 |
| 8 | 11380 |

Ti

pack 3pcs



Ideal for the fixing of meshes and plates.

Microscrew SQ17 ø 1.7

| L | |
|----|-------|
| 9 | 11385 |
| 11 | 11386 |
| 13 | 11387 |



pack 3pcs


Emergency microscrew SQ17 ø 2.0

| L | |
|----|-------|
| 11 | 11390 |
| 13 | 11391 |



pack 3pcs

Ideal for the fixing of bone grafts.



Biomaterials

Granules Adbone BCP



| dimension | g | |
|--------------------|----------|------------|
| Fine 0,1-0,5 mm | 0.5 | BCP010505G |
| | 0.5-5pcs | BCP010505P |
| | 1 | BCP010510G |
| | 1-5pcs | BCP010510P |
| Medium 0,5-1 mm | 0.5 | BCP050105G |
| | 0.5-5pcs | BCP050105P |
| | 1 | BCP050110G |
| | 1-5pcs | BCP050110P |
| Course 1-2 mm | 1 | BCP010210G |
| | 1-5pcs | BCP010210P |



Adbone BCP is a porous synthetic biomaterial in granules, made up of 25 % of Tricalcium phosphate (TCP) and 75% of Hydroxyapatite (Hap). The biphasic composition makes it possible to achieve optimal reabsorption in two stages, compatible with the rapid formation of the bone and the maintenance of the architecture of soft tissues.

Membrane Tisseos

| dimension | |
|-----------|--------|
| 15x20 | TO1520 |
| 15x25 | TO1525 |
| 20x30 | TO2030 |
| 30x40 | TO3040 |

Tisseos is a double layer synthetic membrane, biocompatible and totally reabsorbable, ideal for guided regeneration of the bone and soft tissues.





► Warnings and sales conditions

Warnings

1. Manufacturer responsibility (according to the 93/42 EEC Directive and subsequent amendments)

The Way implant-restoration system is made up of a number of medical devices for Dentistry according to the Directive, aimed at dental restoration of the oral cavity of human beings. The instruments and components dedicated for this purpose make up an integral and indispensable part of the system and must therefore always be used for the application of Way dental implants, scrupulously following the instructions and recommendations supplied by the manufacturer (according to the Directive). Every use of the Way system which is different from the one stated or the use of instruments or components in a manner different to the one foreseen or the use of instruments or components which do not belong to the system, produced by third parties, compromises the functionality of the Way system and is considered as 'improper use', exonerating the manufacturer from any obligation or responsibility. Information concerning the use of Geass products is supplied to the user in written form in paper documentation, like the instructions for use, surgical and restorative protocols, in electronic form (audiovisual and IT instruments) or potentially through practical demonstrations (training courses). These correspond to the current state of art recognized on commercialization of the product and only constitute a supplement to a professional education and experience, as they are not sufficient for an immediate use of the Geass implant systems.

2. User responsibility

Choice and application of the product are acts carried out by the Clinician in total autonomy of judgement and according to the knowledge assumed by the acceptance into the medical-health profession and subsequent professional refreshers; no responsibility can be attributed to Geass for damages of a nature that derive from such acts. The availability of technical-scientific information supporting the client, in fact, does not exonerate the user from the obligation to personally verify the suitability of the product to the purpose of the foreseen procedures. The user is obliged to continually update his knowledge on the development and the applications of the Geass implantological systems. Any use of the system different from the one given, is considered as 'improper use', exonerating the manufacturer from any obligation or responsibility. For uses not expressly foreseen or advised, the user must contact the manufacturer and obtain explicit authorization. The working, handling, and application of the product is performed outside of the manufacturer's control and therefore the responsibility falls to the user. For endoral application of medical devices, it is advisable to always adopt the necessary precautions (e.g. dental dam) in order to eliminate the risk of accidental inhalation.

3. Guarantee

The manufacturer, within the terms and conditions of sale, guarantees that the

products do not have any defects. Geass recognizes a guarantee of twelve months from the delivery date of the product. Geass is obliged to substitute the quantity of products recognized as defective due to manufacture or origin. The guarantee is forfeit and any form of recompense from the manufacturer is excluded should there be improper product use, according to the cases listed in paragraph 1 (manufacturer responsibility) and 2 (user responsibility). Returns must be previously agreed on with the manufacturer and accompanied by the specific documentation. Information on the existence of patents, brand protection rights or other intangible goods is not legally binding.

4. Documentation

The brochures and detailed instructions for use for the implantological Geass systems must be requested from our commercial representatives, area dealers or directly from the head office. Customer service: telephone: +39 0432 669191 – fax +39 0432 665323 e-mail: servizioclienti@geass.it website: www.geass.it Information herein contained shows the state of the art at the moment of commercialization of the product. This does not exonerate the user from the responsibility of personally verifying that the product is suitable for the purposes and procedures foreseen.

5. Seminars and educational course

Geass regularly organizes seminars and educational courses in order to allow users of their products to be informed and refresh their knowledge on the characteristics and on the suitable use of the Geass implant systems.

6. Product identification

All Geass products are identifiable by the article and lot code shown on the accompanying label of the medical devices..

7. Sales packaging

Unless otherwise indicated in the catalog, each product unit identified by the article code is sold in single packaging.

8. Delivery and availability

Geass products are sold to Dentists and Dental laboratories, or for them, according to the relevant competences. Some components may not be available in some Countries or commercial areas.

9. Copyright

Way is a registered brand.

10. Note

For anything not shown in these warnings, see the technical specifications, conditions of use and instructions contained in the Geass informative materials.

Ordering method

1. On placing orders, always refer to the article code.

2. Orders that are received before 12.30 p.m. will be delivered by the end of the following day depending on entity, availability and particular zones. Sales conditions

1. These terms and conditions of sale are intended as accepted by the client on delivery of the order. Any variations, the stipulation of which are hereby illustrated, shall only be valid if accepted by Geass in writing.

2. Regarding market conditions, Geass

reserves the right to modify products, contents of catalogs and prices at any time and with no prior forewarning.

3. Freight charges are paid by the customer. Goods are shipped at the customer's risk even when delivered DAP destination.

4. The delivery terms may undergo variations. Any misunderstandings owing to shipping inefficiency cannot be attributed to Geass.

5. Geass reserves the right to carry out partial delivery.

6. The price list applied is the one valid at the time of the order. Payment of orders must be according to payment method and within the terms established. In the case of default, Geass reserves the right to vary the conditions of payment for subsequent supplies or to put into practice every effective or precautionary measure to totally recoup any outstanding credit.

7. Any complaints, relative to a lack of adherence to the terms and conditions of sale, must be communicated in writing to Geass Customer Service within 8 (eight) days of receiving the goods.

8. Geass srl offers you the possibility to substitute products purchased under the following conditions:

- product cost equal to or above (payment of any difference by client);
- within 12 months of the invoice date and within 6 months of the product going out of date

- date shown on label;

- residual product whole; original packaging complete and sealed;

- product accompanied by transport documentation and a copy of the purchase invoice;

- should these above mentioned conditions not be fulfilled, the product will not be considered suitable and will be returned to the sender and all shipping costs will be charged. Geass srl recognizes the right of withdrawal within 14 working days from the date of the delivery of the goods.

9. Geass declines any responsibility for any involuntary errors in the catalog and price lists.

10. For anything not expressly foreseen in the general terms and conditions of sale, Italian law will be applied. For any disputes, the Court of Udine (Italy) is the competent body.

Document validity

This document substitutes the previous edition.

It is absolutely forbidden to reproduce, even partially, these materials (text and illustrations) without written authorization from Geass S.r.l.



100% producers of smile

Geass is the Italian company which has achieved quality and **innovation** over the last thirty years, offering implant-restorative solutions in order to obtain excellent results.

An internal production line, rigid quality control and next generation technology are all guarantees for reliability of the product, safety and innovation.

www.geass.it

Keep updated with all the latest news on Geass website.



shop.geass.it

A modern tool for your purchases and lots of new promotions.



Social media

News and updates on Geass products and events.



Communication to patient

Materials to help clinicians explain implant treatments and the digital technologies used in clinic.



Smile club

Smile club is a loyalty program that offers access to exclusive benefits, giving even more value to clinicians and their activity.



Geass srl

Via Madonna della Salute, 23
33050 Pozzuolo del Friuli (UD) I
tel+390432669191 info@geass.it
www.geass.it

