

# Product Catalog.

Dental resins, equipment & accessories.



# TCS

An ISO 13485 Certified Company

TCS is a dental manufacturing company dedicated to providing dental laboratories with exceptional quality materials, equipment, and hands-on training courses that are accessible and reliable. At the core of our brand, we believe that one of the key components to creating successful removable appliances is through forming strong partnerships, because for TCS, manufacturing is more than just a business, it is a collectivist culture.

2000

FOUNDED IN CA, USA

3

GENERATIONS

70+

COUNTRIES WE SELL TO



## The TCS Experience

### Our partnership at a glance

**Dynamism.** An energetic and passionate team moves forward together and constantly progresses to add value to our customers.

**Education.** Outstanding technical support and hands-on training designed to grow confidence in creating exceptional TCS dental appliances.

**Innovation.** A philosophy to continuous improvement guides us to develop better methods, tools, materials, educational programs, etc.

**Manufacturing In-House.** The majority of our products, including our resins and equipment are manufactured in-house in CA, USA.

**Together We Smile - TCS Charitable Events.** Our industry is fortunate to be part of a community where giving back literally means putting smiles on people's faces.

**Win - Win Philosophy.** Your success is our success. We aim for long term relationships based on mutual trust, respect and benefit.



# Table of Contents

Unbreakable™ .....	6
iFlex™ .....	10
Karadent™ .....	14
Digital Discs.....	18
Injection Units.....	20
Basic Kits.....	21
Complete Kits.....	22
JP90.....	26
Digital Furnace.....	28
Duplicating Unit.....	29
Waxes / Flasks.....	30

Spatulas / Knives.....	31
Processing Tools.....	32
Finishing Tools.....	33
Polishing Tools.....	34
TCS Fresh.....	36
Get More With TCS.....	38
Technical Support.....	40
Rental Equipment.....	41
Hands-On Training.....	42
Promote Your Department.....	44
Together We Smile.....	46

# Unbreakable™


Flexible Partial



**Monomer Free**  
**Hypoallergenic**



- Extremely resilient, high performance nylon thermoplastic.
- Guaranteed unbreakable due to its unique properties.
- Designed for exceptional memory to hold its form.
- Offers an ideal degree of flexibility.
- A modern fit, lightweight and durable.




# Unbreakable™

TCS resins are available in bulk or aluminum cartridges that come in various sizes/diameters.

**Be sure to order the correct diameter cartridge to fit your injection unit.**

*Refer to [youtube.com/tcsdentalinc](https://www.youtube.com/tcsdentalinc) (TCS Short Videos Playlist).*

**Millable discs available for CAD/CAM - Refer to pages 18-19.**



Shade	Size	25 mm	28 mm	22 mm
Light Pink (Shade 1)	S	4002-01/S1		
	M	4000-01/M1	4000-01/MF1	4005-01/M22-1
	L	4005-01/L1	4005-01/LF1	4005-01/L22-1
Standard Pink (Shade 2)	S	4002-01/S2		
	M	4000-01/M2	4000-01/MF2	4005-01/M22-2
	L	4005-01/L2	4005-01/LF2	4005-01/L22-2
Light / Dark Pink (Shade 3)	S	4002-01/S3		
	M	4000-01/M3	4000-01/MF3	4005-01/M22-3
	L	4005-01/L3	4005-01/LF3	4005-01/L22-3
Dark Pink (Shade 4)	S	4002-01/S4		
	M	4000-01/M4	4000-01/MF4	4005-01/M22-4
	L	4005-01/L4	4005-01/LF4	4005-01/L22-4
Natural (Shade N)	S	4002-01/SN		
	M	4000-01/MN	4000-01/MFN	4005-01/M22-N
	L	4005-01/LN	4005-01/LFN	4005-01/L22-N

**Each pack includes the following:**

5 delivery bags | 5 patient care instructions | 5 samples of TCS Dental Appliance Cleaner (1-week supply).\*

**Each kilo includes the following:**

50 delivery bags, 50 patient care instructions.

*Not applicable to JP Bullets. USA only.\**




25 mm w / tab	25 mm w / nozzle	1 Kilo	500 Grm	JP Bullets
				4002-JP/S1
4005-01/MT25-1		4003-01B/1	4003-01B/1-.5	
4005-01/LT25-1	4005-01/LNZ25-1			4005-JP/L1
				4002-JP/S2
4005-01/MT25-2		4003-01B/2	4003-01B/2-.5	
4005-01/LT25-2	4005-01/LNZ25-2			4005-JP/L2
				4002-JP/S3
4005-01/MT25-3		4003-01B/3	4003-01B/3-5	
4005-01/LT25-3	4005-01/LNZ25-3			4005-JP/L3
				4002-JP/S4
4005-01/MT25-4		4003-01B/4	4003-01B/4-.5	
4005-01/LT25-4	4005-01/LNZ25-4			4005-JP/L4
				4002-JP/SN
4005-01/MT25-N		4003-01/BN	4003-01/BN-.5	
4005-01/LT25-N	4005-01/LNZ25-N			4005-JP/LN

**iFlex™**  
Flexible Partial

**BPA  
FREE**  
**Monomer Free  
Hypoallergenic**



- Highly translucent, polyolefin thermoplastic.
- Simple to finish and polish.
- Stain resistant and easy to maintain.
- Clinically unbreakable.
- Offers an ideal degree of flexibility.
- A modern fit, lightweight and durable.




# iFlex™

TCS resins are available in bulk or aluminum cartridges that come in various sizes/diameters.

**Be sure to order the correct diameter cartridge to fit your injection unit.**

*Refer to [youtube.com/tcsdentalinc](https://www.youtube.com/tcsdentalinc) (TCS Short Videos Playlist).*

*Millable discs available for CAD/CAM - Refer to pages 18-19.*



Shade	Size	25 mm	28 mm	22 mm
Light Pink (Shade 1)	S	5002-01/S1		
	M	5000-01/M1	5000-01/MF1	5005-01/M22-1
	L	5005-01/L1	5005-01/LF1	5005-01/L22-1
Standard Pink (Shade 2)	S	5002-01/S2		
	M	5000-01/M2	5000-01/MF2	5005-01/M22-2
	L	5005-01/L2	5005-01/LF2	5005-01/L22-2
Light / Dark Pink (Shade 3)	S	5002-01/S3		
	M	5000-01/M3	5000-01/MF3	5005-01/M22-3
	L	5005-01/L3	5005-01/LF3	5005-01/L22-3
Dark Pink (Shade 4)	S	5002-01/S4		
	M	5000-01/M4	5000-01/MF4	5005-01/M22-4
	L	5005-01/L4	5005-01/LF4	5005-01/L22-4
Natural (Shade N)	S	5002-01/SN		
	M	5000-01/MN	5000-01/MFN	5005-01/M22-N
	L	5005-01/LN	5005-01/LFN	5005-01/L22-N

**Each pack includes the following:**

5 delivery bags | 5 patient care instructions | 5 samples of TCS Dental Appliance Cleaner (1-week supply).\*


**Each kilo includes the following:**

50 delivery bags, 50 patient care instructions.

*Not applicable to JP Bullets. USA only.\**



25 mm w / tab	25 mm w / nozzle	1 Kilo	500 Grm	JP Bullets
--				5005-JP/S1
5005-01/MT25-1		5003-01B/1	5003-01B/1-.5	--
5005-01/LT25-1	5005-01/LNZ25-1			5005-JP/L1
--				5005-JP/S2
5005-01/MT25-2		5003-01B/2	5003-01B/2-.5	--
5005-01/LT25-2	5005-01/LNZ25-2			5005-JP/L2
--				5005-JP/S3
5005-01/MT25-3		5003-01B/3	5003-01B/3-.5	--
5005-01/LT25-3	5005-01/LNZ25-3			5005-JP/L3
--				5005-JP/S4
5005-01/MT25-4		5003-01B/4	5003-01B/4-.5	--
5005-01/LT25-4	5005-01/LNZ25-4			5005-JP/L4
--				5005-JP/SN
5005-01/MT25-N		5003-01B/N	5003-01/BN-.5	--
5005-01/LT25-N	5005-01/LNZ25-N			5005-JP/LN



# Karadent™

Full Dentures  
Clear Frameworks



**Monomer Free**  
**Hypoallergenic**



- Semi-rigid, microcrystalline polymer thermoplastic.
- Ideal for full dentures, semi-flexible partial dentures, clear clasps, frameworks, and nightguards.
- An upgrade from traditional acrylic dentures yet as simple to finish and polish.
- Finishes completely transparent in our Clear shade.
- A modern fit, lightweight and durable.






# Karadent™

TCS resins are available in bulk or aluminum cartridges that come in various sizes/diameters.

**Be sure to order the correct diameter cartridge to fit your injection unit.**

*Refer to [youtube.com/tcsdentalinc](https://www.youtube.com/tcsdentalinc) (TCS Short Videos Playlist).*

**Millable discs available for CAD/CAM - Refer to pages 18-19.**



Shade	Size	25 mm	28 mm	22 mm
Standard Pink (Shade 2)	S	2000-01/S2	--	--
	M	2000-01/M2	2000-01/MF2	2005-01/M22-2
	L	2005-01/L2	2005-01/LF2	2005-01/L22-2
Dark Pink (Shade 4)	S	2000-01/S4	--	--
	M	2000-01/M4	2000-01/MF4	2005-01/M22-4
	L	2005-01/L4	2005-01/LF4	2005-01/L22-4
Natural (clear) (Shade N)	S	2000-01/SN	--	--
	M	2000-01/MN	2000-01/MFN	2005-01/M22-N
	L	2005-01/LN	2005-01/LFN	2005-01/L22-N

**Each pack includes the following:**

5 delivery bags | 5 patient care instructions | 5 samples of TCS Dental Appliance Cleaner (1-week supply).\*


**Each kilo includes the following:**

50 delivery bags, 50 patient care instructions

*Not applicable to JP Bullets. USA only.\**



25 mm w / tab	25 mm w / nozzle	1 Kilo	500 Grm.	JP Bullets
--	--			2002-JP/S2
--	--	2003-01/B2	2003-01B/2-.5	--
--	2005-01/LNZ25-2			--
--	--			2002-JP/S4
--	--	2003-01/B4	2003-01B/4-.5	--
--	2005-01/LNZ25-4			--
--	--			2002-JP/SN
--	--	2003-01/BN	2003-01/BN-.5	--
--	2005-01/LNZ25-N			--




## *Dive Into the Digital Age*



***Available Soon***

## ***Millable Discs*** *for a modern fit*



- Digital technology allows for ultimate precision.
- Appliances come out with a smooth texture, close to polishing.
- Polishing process is quick (approx. 2 minutes).
- Milling an appliance from a disc is kinder on your tools.

***Contact us for details***  
**(866) 426-2970**

## TCS Injection Units.

TCS offers a variety of Injection Units to fit laboratories preferences and budget.

**MOST POPULAR**



**Hydraulic Injector (electric).**

- Equal pressure with every injection.
- Minimal maintenance.
- No air connections required.
- Requires an external heater (TCS Digital Furnace)

Injector only 110V	3504-01
Injector only 220V	3504-02

**SPECS.**

Yields approx. 30 cases per day (up to 120 cases with additional furnaces).

**H:** 22" **L:** 11" **W:** 6.5"  
50 lbs. (23 kilos).

**2 year warranty**  
**Made in USA**

**MOST ECONOMICAL**



**Manual Injector.**

- Minimal maintenance.
- Built to last.
- Simple to use.
- Requires an external heater (TCS Digital Furnace)

Injector only 110V	3506-01
--------------------	---------

**SPECS.**

Yields approx. 30 cases per day (up to 120 cases with additional furnaces).

**H:** 22" **L:** 22" **W:** 17"  
42 lbs. (19 kilos).

**2 year warranty**  
**Made in USA**

**ALL IN ONE**



**Automatic Injector.**

- Compact design with a built-in furnace.
- Equal pressure with every injection.
- Injects with a touch of a button.
- Requires hookup to air pressure (min 140 psi)

Injector only 110V	3525-01
Injector only 220V	3525-02

**SPECS.**

Yields approx. 30 cases per day.

**H:** 23" **L:** 8" **W:** 13.5"  
40 lbs. (18 kilos).

**1 year warranty**  
**Made in USA**

**USE 25MM CARTRIDGES FOR ALL INJECTORS**

For TCS Digital Furnace, see page 28


**YouTube:** [youtube.com/tcsdentalinc](https://www.youtube.com/tcsdentalinc)  
(TCS Equipment playlist)

## Basic Kits.

Kits include the basic essentials to inject TCS Resins.

**Hydraulic Basic Kits.**

110V	3511-01
220V	3511-02



**Manual Basic Kits.**

110V	3508-01
220V	3508-02



1 digital furnace

1 cylinder with handle

1 knock out rod  
1 knock out stand

1 thick bronze plug  
1 thin bronze plug

1 injection flask




2 cartridges of each:  
Unbreakable™, iFlex™,  
Karadent™



1 mold release spray

**Automatic Air Basic Kits.**

110V	3528-01
220V	3528-02



1 injection flask



2 cartridges of each:  
Unbreakable™, iFlex™



1 mold release spray


Note: It is not recommended to inject Karadent™ by tcs® using the Automatic Air Injector.

## TCS Complete Injection Kits.

Includes all the essentials necessary to process TCS Resins.

### Hydraulic Complete Kits.

110V 3501-01  
220V 3501-02



### Manual Complete Kits.

110V 3505-01  
220V 3505-02



2 pks - Unbreakable™  
2 pks - iFlex™  
2 pks - Karadent™



1 digital furnace



2 cylinders with handle



1 knock out rod  
1 knock out stand



2 thick bronze plug  
1 thin bronze plug



2 injection flasks



10 pcs. palate wax



10 pcs. lingual wax



10 pcs. sprue wax



3 pcs. 9.5 sprue wax



1 carving knife



1 cement



1 EZ flow sample



1 mold release spray



1 fusing liquid



1 hot air welder




Twist drill samples



TCS putty (400 gms)



1 econo cutter  
1 econo mandrel



1 white knock down  
wheel & chuck



1 rubberized abrasive 304  
1 rubberized abrasive 306



10 pink grinding wheels



10 brown rubber wheels



5 green mounted stones



2 large rubber points  
2 small rubber points



2 mini fuzz buffs



1 pink zirconia coarse bur



1 lb polycril



1 turbo buff



4 rag buffs



sample pieces of  
polishing paste




1 sample model  
upon request

















Items in the kits may change without notice. We strive for continuous improvement, if we discover a product that works better with our resins, we may switch out an item(s).



## TCS Complete Injection Kits.

Includes all the essentials necessary to process TCS Resins.




 3 pks - Unbreakable™ 3 pks - iFlex™	 2 injection flasks	 10 pcs. palate wax	 10 pcs. lingual wax
 10 pcs. sprue wax	 3 pcs. 9.5 sprue wax	 1 carving knife	 1 cement
 1 EZ flow sample	 1 mold release spray	 1 fusing liquid	 1 hot air welder
 Twist drill samples	 TCS putty (400 gms)	 1 econo cutter 1 econo mandrel	 1 white knock down wheel & chuck

 1 rubberized abrasive 304 1 rubberized abrasive 306	 10 pink grinding wheels	 10 brown rubber wheels	 5 green mounted stones
 2 large rubber points 2 small rubber points	 2 mini fuzz buffs	 1 pink zirconia coarse bur	 1 lb polycril
 1 turbo buff	 4 rag buffs	 sample pieces of polishing paste	 1 sample model upon request

*Note: It is not recommended to inject Karadent™ by tcs® using the Automatic Air Injector.*

Items in the kits may change without notice. We strive for continuous improvement, if we discover a product that works better with our resins, we may switch out an item.



## JP90 Handheld Injector.



Used in conjunction with the TCS Digital Furnace. Sold separately.

### Injects:

- » 1-3 tooth unilaterals.
- » Up to 5 consecutive tooth partials/saddles, saddles for combination cases.
- » Repairs.

JP90 Injector only      7000-JP

### SPECS.

H: 6" L: 2-1/4" W: 11"  
2 lbs. (1 kilo).

**1 year warranty**  
Made in USA

▶ **YouTube:** [youtube.com/tcsdentalinc](https://www.youtube.com/tcsdentalinc)  
(TCS Equipment playlist)



### JP Flask.

Fabricated from high precision stainless steel to withstand warpage under extreme heat.

7005-JP



### Metal Sprues.

Used to create a gateway for injection.

7014-JP



### Matrix Pusher.

Used to facilitate putty application.

7010-JP



### TCS Putty.

Used to create a matrix when injecting with the JP90 Handheld Injector. Can also be used to invest for clean injections.

Base & Catalyst      4296-01  
400gms ea.

Base & Catalyst      4291-01  
1 kilo ea.

Base & Catalyst      4292-01  
5 kilo ea.

## JP90 Starter Kits.

**JP90 Starter Kits.**  
JP Starter Kit      7001-JPU  
JP Starter Kit with Digital Furnace      7020-JPU

Includes TCS Digital Furnace



1 JP90 handheld injector

JP bullets (small)  
1 pk - Unbreakable™  
1 pk - iFlex™

JP bullets (large)  
1 pk - Unbreakable™  
1 pk - iFlex™

JP flask

TCS putty  
400 gms


2 waxed cases for practice

1 matrix pusher  
5 metal sprues

## JP90 Bullets.

Refer to the table below for bullet shades, sizes, and item codes.

Unbreakable™	page 9
iFlex™	page 13
Karadent™	page 17



## Digital Furnace & Accessories.



- » Simple and practical.
- » Melts thermoplastic resin.
- » Compact design.
- » Reliable / low maintenance.

110V	3512-01
220V	3512-02

**SPECS.**  
**H:** 6" **L:** 1.5" **W:** 11"  
 1.3 lbs. (0.6 kilos).

**2 year warranty**  
 Made in USA

▶ **YouTube:** [youtube.com/tcsdentalinc](https://youtube.com/tcsdentalinc)  
 (TCS Equipment playlist)



### Melting Cylinder with Handle.

Fabricated from high precision stainless steel to withstand warpage under extreme heat.

**3703-01**



### Knock Out Stand and Rod.

Steel tools used for purging the melting cylinder from the remaining cartridge and bronze plug after injection.

**Knock Out Stand. 3705-01**

**Knock Out Rod. 3707-01**



### Bronze Plugs.

Thick bronze plug is used when melting medium and small cartridges in the furnace.  
 Thin bronze plug is used when melting large cartridges in the furnace.

**Bronze Plug Thin. 3709-01**

**Bronze Plug Thick. 3712-01**

## Duplicating Unit - 2 Gallon. Model TCS105



- » Double boil-out tank helps prolong the life of hydrocolloid paste.
- » No water connections required.
- » Simple design and easy to use.

110V	3614-01
Replacing previous unit 3611-01	

**SPECS.**  
 Yields approx. 18 pour-ups (depending on flask size).  
**H:** 21" **L:** 14-1/4" **W:** 10"  
 20 lbs (9 kilos).

**1 year warranty**  
 Made in USA



**No frills unit ...**  
*For value and reliability*

## Waxes / Flasks.

### Palate Wax.

Save time waxing with pre-formed palates created with the recommended thickness of a TCS palate.

5015-01 80/box



### Lingual Wax.

Pre-formed lingual bars are ideal for clasps, flange, and lingual connector areas.

5020-01 1 lb



### Sprue Wax.

Easy to work with, does not get brittle. Ideal diameter for spruing. 5 gauge.

5025-01 1 lb




### Sprue Wax 9.5.

Connect your sprue wax to the 9.5 Sprue Wax only at the tip of the flask opening to avoid unwanted plaster sticking to it.

5027-01 5 pcs

5026-01 1 lb (~30 pcs)

Each piece yields ~10 injections.



### Injection Flask.

Milled from one piece of aluminum block. Extra smooth finish facilitates the de-flasking process. Includes a special knock-out button to simplify case removal and prolong the life of the flask.

3701-01



### Pour Flask.

Used for acrylic pour technique and/or to duplicate models. Manufactured in light-weight aluminum. Easy to assemble, disassemble and clean.

3650-01



### Nightguard Pour Flask.

Used for the fabrication of nightguards and/or duplicate models. Manufactured in light-weight aluminum. Easy to assemble, disassemble and clean.

3651-01



### Brass Spruing Tools.

Used to sprue hydrocolloid paste for acrylic pour technique.

Medium (12mm dia) 3653M-01

Small (8.75mm dia) 3653S-01



## Spatulas / Knives.

### Plaster Spatula.

German steel blade with wooden handle. 4-1/4" flexible blade.

4230-01



### Plaster Knife.

German steel blade with wooden handle. 2-1/4" blade.

4231-01



### Fahnstock Knife and Spatula.

German steel blades with wooden handle for waxing.

4232-01



### Wax Spatula.

German steel double ended spatula for waxing.

4233-01



### Wax Spatula.

German steel double ended spatula for waxing.

4234-01



### Wax Carver.

German steel double ended spatula and carver for waxing.

4235-01



### TCS Carving Knife.

Excellent for removing extended fibers that sometimes appear when processing TCS partials. Also used to define the necklines of the TCS restorations when necessary.

4201-01



### Replacement Blades.

Replacement blades for TCS Carving Knife.

4202-01 12 pcs





## Processing Tools.

### TCS Cement.

To secure teeth in the stone prior to injection. Mild enough to keep teeth from bonding to the stone during injection.

3801-01 1 oz



### Fusing Liquid.

For repairs and additions to existing TCS restorations. The Fusing Liquid is ONLY used when flasking and injecting are necessary.

3812-01 1/2 oz



### Mold Release Spray.

Colorless and non-staining multi-purpose lubricating spray, used to ease the purging of the furnace accessories and TCS cartridges during injection process.

3814-01 12 oz



### Hot Air Welder.

Hand torch in a slim design for minor repairs or additions to TCS restorations. Refillable with butane. (Butane not included).

3714-01



### Twist Drills.

To create the diatoric retention in acrylic teeth necessary for fabricating TCS restorations.

005 5005-01 6 pcs

006 5006-01 6 pcs

009 5009-01 6 pcs

012 5012-01 6 pcs



### EZ Flow.

Separating medium for processing thermoplastics. Eliminates stone sticking to the partial, leaving it shiny and clean.

4217-01 16 oz

4218-01 1 gallon



*\* Not recommended when hot flask injection is required, e.g. Karadent™.*

### Empty Cartridges.

Empty cartridges are sold in packs of 50.

	25mm	28mm	22mm
Small	3003-01S		
Medium	3003-01M	3003-01MF	3003-01M22
Large	3003-01L	3003-01LF	3003-01L22



## Finishing Tools.

### Econo Cutter XL.

Use to cut off sprues from TCS restorations after injection.

4203-01



### Econo Cutter Mandrel.

Mandrel required for the econo cutter.


4204-01



### White Knock Down Wheel.

One of the first steps in the finishing stage, cuts excess bulk and trims partial to the desired length.

4211-01



### White Knock Down Chuck.

Mandrel used with the White Knock Down Wheel and Rubberized Abrasives.

4212-01



### Rubberized Abrasives.

Long-lasting wheel for smoothing the partial during the finishing stage.

Thin 4304-01

Thick 4306-01



### Pink Grinding Wheels.

Excellent and economical product for the finishing stage.

4205-01 100/box



### Brown Rubber Wheels.

Use prior to pumicing in order to obtain a high shine. It is crucial that the entire partial is rubberized and smoothed.

4503-01 100/box



### Mini Fuzz Buff.

Use after the Brown Rubber Wheel to further smoothen partial prior to pumicing.

4222-01 12 pcs



## Finishing Tools / Polishing Tools.

### Brown Rubber Points.

Ideal to eliminate scratches and smoothen the surface.

Small	4513-01	12 pcs
Large	4514-01	12 pcs



### Green Mounted Stones.

Used during the finishing stage. Cuts thermoplastics moderately, and leaves a smooth surface. Ideal for minor chairside adjustments.

4501-01	100/box
---------	---------



### Carbide Bur.

Specifically designed to cut thermoplastics. Long-lasting with high cutting efficiency.

50G79-01

50251G-01



### Pink Zirconia Fine Bur.

Use during the finishing stage to leave the surface smooth and prepped for the Brown Rubber Wheel. Also, ideal for chairside adjustments. Max RPM=12,000. Use with a hand piece.

5035Z



### Pink Zirconia Coarse Bur.

Use during the finishing stage to cut excess material and trim appliance to the desired length. Gross reduction. Max RPM=16,000. Use with a hand piece.

5045Z



### Polycril.

An advanced formula for pumicing acrylics and thermoplastics at high-speed.


4298-01	1 lb
---------	------



### Turbo Buff.

Pumicing buff with siliconized cotton cloth. Use for pumicing only.

4208-01



### Rag Buff.

Stitched only in the center for the extra softness needed for polishing TCS restorations.

4207-01



## Polishing Tools.

### Mini Rag Buff.

Recommended to dentists for polishing touch-ups after minor adjustments on TCS appliances. Use with Green Hi-Shine polishing paste.

4221-01	12 pcs
---------	--------



### Brown Tripoli.

Use after pumicing to prepare appliance for polishing. Effective for all TCS resins. Use with a dry Rag Buff.

4206-01	2 lbs
---------	-------



### iShine.

Use after Brown Tripoli to polish surface of iFlex™ partials (and other polyolefin thermoplastics). Use with a dry Rag Buff.

4220-01	10.75 oz
---------	----------



### Green Hi-Shine.

Use after Brown Tripoli to polish surface of Karadent™ and Unbreakable™ partials. Use with a dry Rag Buff.

4209-02	1 lb
---------	------



### Sample Model Unbreakable™.

An ideal marketing tool for your laboratory.

Upper Partial	4011-01U
---------------	----------

Lower Partial	4011-01L
---------------	----------



### Sample Model: Karadent™.

An ideal marketing tool for your laboratory.

Upper Full Denture	4016-01U
--------------------	----------

Lower Partial Denture	4017-01L
-----------------------	----------

Upper Clear Framework	4018-01U
-----------------------	----------



### Sample Model: iFlex™.

An ideal marketing tool for your laboratory.

Upper Partial	4009-01U
---------------	----------

Lower Partial	4009-01L
---------------	----------



### Dr. Adjustment Kit.\*

Ideal for quick and easy chairside adjustments.

4014-02



\*Includes: 1 Green Point, 1 large Rubber Point, 1 small Rubber Point, 1 Pink Zirconia Fine Bur, 1 Carving Knife, 1 sample of Green Hi-Shine, 2 Mini Rag Buffs, 5 Mini Fuzz Buffs.



**TCS includes patient care instructions and a sample pack of denture cleaner with each cartridge.\***

Include them with the finished case so that the patient can maintain a clean and vibrant appliance.

\*USA only



A clean appliance is vital to patients' health and the appearance and longevity of dental appliances. TCS Concentrated Dental Appliance Cleanser is a safe solution for quick and effective results.



**TCS Concentrated Dental Appliance Cleaner.**

Professional strength cleanser. Deep cleans, removes stains and odors. Ideal for most types of removable dental appliances.

4505-02      24 pcs  
6 month supply

**TCS Dental Appliance Cleaner.**

TCS Dental Appliance Cleaner sample packets. Box of 75 packets (1-week supply ea.).

4505FR-SMP      75 pcs

**TCS Dental Appliance Brush.**

Resilient soft bristle brush to thoroughly clean dental appliances without scratching the surface.

4506FR-01

FREE DISPLAY WITH ORDER OF 9+ BOXES



# Get More with TCS

- Technical Support.
- Equipment Rental.
- Training.
- Promotional Material.
- TCS Gives Back.

## *We Value Your Success*

Quality tools and support help improve your laboratory's productivity and quality when fabricating a TCS Dental Appliance.

**Connect with us to stay up to date on our latest information, videos and techniques.**

**Subscribe to our newsletter:**

[www.tcsdentalinc.com](http://www.tcsdentalinc.com)

**Follow Us:**



@tcsdentalinc



# Technical Support.

How to find the most recent videos on our You Tube channel and register for TEXT Support.


## YouTube Channel: *tcsdentalinc*

Subscribe to view short clips that will help with production and quality.

### Playlists:

#### QUICK TIPS

1-3 min videos



#### DESIGN LIBRARY

Design & blockout



#### SHORT VIDEOS

TCS Staff shares news



#### EQUIPMENT

Animated equipment demonstration




## Text Support.

Register your phone for quick text support.

Text your name, city & state to (562) 241-1225.

## Phone.

(866) 426-2970 | (562) 426-2970



TECHNICAL QUESTIONS?

Text us!  
(562) 241-1225

# Rental Equipment Program.

## On the fence?

Rent a TCS injection unit for 20 days prior to investing.



Try Before You Buy.

Contact us to get started:

sales@tcsdentalinc.com

## Hands-On Training.

### ***We are dedicated to continuing education.***

In our comprehensive hands-on training classes technicians will learn every step needed to fabricate a flexible partial. Courses cover the importance of case design, quick finishing, high shine polishing and everything in between.

***We highly recommend that new and experienced technicians attend our hands-on training classes to stay up to date on the latest technical practices.***

**For an Up-To-Date Schedule & Details.**

**Visit:** [www.tcsdentalinc.com/pages/training](http://www.tcsdentalinc.com/pages/training)

**Contact us:** [sales@tcsdentalinc.com](mailto:sales@tcsdentalinc.com) | [www.tcsdentalinc.com](http://www.tcsdentalinc.com)

Contact us for  
information on

SEMINARS

LUNCH N' LEARNS

WEBINARS

TRAINING VIDEOS

TRAINING MODULES

REQUESTING A VIDEO

*"The presentation and lecture was the best I've ever attended. Well worth taking the day for continuing education. I thoroughly enjoyed the day." ~ TX Training*

*"I learned a lot of new things, including tips on proper case design, time saving techniques, and troubleshooting." ~ CA Training*

*"Training course was off-the-charts." ~ NY Training*



# Promote Your Removable Department.

Grow your removable department by utilizing these simple marketing tools.



## Digital Media.

In today's digital age one of the most effective ways to promote your removable department is through online content. Enhance your website / social media platforms by embedding and sharing our digital marketing materials.

Contact us or visit our website for updated videos.

Ask about our banner / social media look book to utilize digital ads and more.

# Order Free Promotional Material



## Patient Education Brochures.

Informative brochure for patients to learn about TCS appliances as a solution for missing teeth. Distribute to your accounts to display in their waiting area.

Size: 8.5" x 11".  
Folded: 8.5" x 3.6".

A-8010-A1



## Dr. Informational Booklet.

Includes everything a doctor needs to know about TCS appliances. E.g. impression taking, insertion, chairside adjustment techniques. A must have for dentists.

Size: 6" x 7.9".  
(19 page booklet).

A021413



## Dr. Flyer.

Convenient one page informative flyer to distribute to your accounts. Includes information on all TCS resins & tools for their practice.

Size: 8.5" x 11".

122010



## Statement Stuffers.

Include these product stuffers in account statements. Includes space for you to add your logo and contact info.

Size: 8.5" x 3.5".

iFlex™ A020613  
Karadent™ A040511  
Unbreakable™ A033116



## Shade Guide.

Provide a shade guide to your accounts as a reference to match patient's natural gingiva to TCS resin shades.

Size: 3" x 2.5".

8015-02

ASK

about our affiliate programs

## JOIN • PROMOTE • EARN

Contact us today for more information.

Pascui Krajewski  
pascui@tcsdentalinc.com

Feel free to contact us for promotional materials, digital material for websites / social media, and assistance with clinical studies.

## CONTACT US TODAY.

Phone: 866-426-2970  
Online: www.tcsdentalinc.com (labs & technicians - order promotional material).

Ships to USA only.

# together we smile

TCS GIVES BACK



Dominican Republic, 2019

TCS travels to remote areas across the globe to deliver partials and dentures to people in need. Our outreach program is possible thanks to our collaboration with participating local dentists and technicians. We are always looking for new opportunities.

**Contact us for possible collaborations:**

[info@tcsdentalinc.com](mailto:info@tcsdentalinc.com)







2619 Lime Ave. Signal Hill, CA 90755, USA

(866) 426-2970 | (562) 426-2970

@**tcsdentalinc**

www.**tcsdentalinc**.com