



DMP DENTAL

RESTORATIVE SOLUTIONS

Easily & confidently create long lasting, natural-looking restorations
that your patients won't stop smiling about!



CLASS I-IV

NANOCERAM-BRIGHT

Nano-hybrid universal composite with a special blend of micro and nano fillers for class I-IV restorations.



CLASS I-IV

NANOCERAM-BRIGHT Flow

Flowable version of NANOCERAM-BRIGHT that is highly thixotropic for easy handling and precise application.



CLASS I-II

BRIGHT Posterior

Designed for posterior class I & II restorations its advanced technology results in remarkable physical properties.

Superior Esthetic Results. Exceptional Physical Properties

NANOCERAM-BRIGHT and NANOCERAM-BRIGHT Flow

With 30+ years of experience; DMP Dental has developed an innovative manufacturing process that integrates a unique blend of micro and nano-fillers into the resin matrix creating an easy to sculpt, nano-hybrid composite that achieves superior physical and aesthetic results.



WHY CHOOSE NANOCERAM-BRIGHT

Superior Aesthetics

- Chameleon effect for natural-looking restoration.
- Average partical size of less than 700nm for high polishability & polish retention
- High color retention due to low water absorption & solubility

Exceptional Handling Characteristics

- Excellent sculptability
- No slumping; stays in place
- Non-sticky and easy to polish.
- Highly thixotropic flowable available

Long-Lasting Results

- Low water absorption & solubility ensuring durability.
- Low polymerization shrinkage resulting in less stress on the bond & better marginal integrity.
- Excellent flexural strength reduces deformation under load-chewing.

BRIGHT Posterior

Composite designed for posterior class I and II restorations. Outstanding handling characteristics ensure fast & easy placement with excellent aesthetic results. Its advanced technology results in remarkable physical properties which enhances its



WHY CHOOSE BRIGHT POSTERIOR

Excellent Strength

- Superior compression strength to withstand chewing forces.
- Lowest polymerization shrinkage resulting in less stress on the bond.
- Fast curing up to 4 mm increments

Exceptional Handling Characteristic

- Creamy consistency allows you to bulk fill the posterior without overfilling.
- Thixotropic flow virtually eliminates the need for heat or vibrating devices.
- Non-sticky
- Strongly radiopaque

Superior Wear Resistance

- Low water absorption & solubility ensuring long-lasting durability and color retention
- High fracture and abrasion resistance reducing cracks or breaks

NANOCERAM-BRIGHT Clinical Case

Replacing a discolored IV restoration with NANOCERAM-BRIGHT performed by Dr. Ross Nash, DDS

Before



After



Step by Step Process



Product Number and Description

NANOCERAM-BRIGHT

REFILL 4g x 2

160111401	shade A1
160111402	shade A2
160111403	shade A3
160111405	shade A3.5
160111406	shade A4
160111407	shade B1
160111408	shade B2
160111409	shade B3
160111411	shade C1
160111412	shade C2
160111413	shade C3
160111416	shade D3
160111404	shade OA2
160111419	shade incisal
160111424	shade bleach extra
160111425	shade bleach

NANOCERAM-BRIGHT

Unit Dose 0.25g x 20

160122501	shade A1
160122502	shade A2
160122503	shade A3
160122505	shade A3.5
160122506	shade A4
160122507	shade B1
160122508	shade B2
160122509	shade B3
160122511	shade C1
160122512	shade C2
160122513	shade C3
160122516	shade D3
160122504	shade OA2
160122519	shade incisal
160122524	shade bleach extra
160122525	shade bleach

NANOCERAM-BRIGHT Flow

2 Syringe Kit 2g x 2, syringe tips

170112301	shade A1
170112302	shade A2
170112303	shade A3
170112305	shade A3.5
170112307	shade B1
170112308	shade B2
170112311	shade C1
170112315	shade D2
170112304	shade OA2
170112319	shade incisal
170112325	shade bleach

BRIGHT Posterior Clinical Case

Class V restoration by Dr. Ross Nash, DDS


Before



After



Step by Step Process



Product Number and Description

BRIGHT Posterior

REFILL 4g

Unit Dose 0.25g x 20







180112501 shade A1
180112502 shade A2
180112503 shade A3

180122501 shade A1
180122502 shade A2
180122503 shade A3

Exceptional physical properties proven

Higher Resistances to Masticatory Forces



NANOCERAM-BRIGHT and BRIGHT Posterior outperform the competition with higher compressive, dimetral, and flexural strengths resulting in a composite that better resists wear and deformation under load-chewing forces. These tests show how well composites withstand forces applied to the top and bottom (Compressive Strength), sides (Diametral strength) and 3-point bend technique (Flexural Strength).



to outperform the competition.

Minimal Shrinkage

Unique blend of micro and nano fillers result in minimal shrinkage. This means there is less stress on the adhesive bond area, lower deformation of the tooth structure and better margin quality.



Superior Esthetic Results proven to outperform the competition.

Chameleon Effect

NANOCERAM-BRIGHT and BRIGHT Posterior mimic the natural tooth with matched refractive indexes for both the fillers and monomers and lifelike translucency which allow our composite to blend into the natural tooth structure and surrounding area giving them a true chameleon effect.





Eight direct veneers #5-#12 using one shade of NANOCERAM-BRIGHT performed by Dr. Ross Nash, DDS



BRIGHT Posterior restoration blends into the surrounding teeth performed by Dr. Ross Nash, DDS


Lowest Water Sorption to Help Maintain Color

NANOCERAM-BRIGHT and BRIGHT Posterior have one of the lowest water sorption rates on the market ensuring that the original color is maintained longer throughout the life of the restoration.



Unsurpassed Polish Retention

NANOCERAM-BRIGHT has an average particle size of less than 700nm for superior polishability and polish retention leading to long lasting natural looking aesthetic results.



What the Experts Say

"(NANOCERAM-BRIGHT) Color matching was excellent and after finishing the material seemed to have a chameleon effect as it blended seamlessly into the surrounding teeth." Dr. Ross Nash, DDS



"NANOCERAM-BRIGHT's chameleon-like properties, along with comprehensive range of restorative tints, offer excellent shade matching for restorations that literally disappear" Dr. George Freedmand, DDS



"I have been using DMP composites for about 2 years full time in my practice. I really like using BRIGHT Posterior because it condenses well, resulting in broad, strong interproximal contacts." Dr. William Jones, DDS



"I found NANOCERAMBRIGHT to be a proven, cost effective nano-hybrid composite that excellent results in any application." Dr. Joseph Pelerin



"When using BRIGHT POSTERIOR I was able to carve most of the anatomy with instruments and saved significant chair time with polishing." Dr. Ross W. Nash, DDS

Product	Composite Type	Class I	Class II	Class III	Class IV	Class V	Sealant	Liner
Nanoceram Bright	Universal Composite	●	●	●	●	●		
Nanoceram Bright Flow	Flowable Composite			●	●	●	●	●
Bright Posterior	Posterior Composite	●	●					



DMP DENTAL

www.dmpdental.com | info@dmpdental.com