

## C339

Rotating anode X-ray mammography tube unit. Available in 3 version with different HV and LV receptacles positions.

## **Specifications**

Overall length Maximum diameter Tube assembly net weight Nominal X-ray tube assembly voltage Maximum tube assembly heat content Maximum continuous heat dissipation Minimum tube assembly inherent filtration Maximum leakage radiation at 1m from focal spots 45 µmGy/h

340 mm 140 mm 13 kg 49 kv 375 kJ (500 kHU) 80 W (6400 HU/min) 0.5 mm-Be

## Insert for Reloading

- XM12
- XM12T
- XM15 • XM15T
- XM1016
- XM1016T

## **3 Versions Available**

- C339E lateral HV and LV sockets, VARIAN version
- C339C coaxial HV and LV sockets, COMET version
- C339V coaxial HV and LV sockets, VARIAN version

