

VOBIS

VOUDOURIS – OPEN – BITE – INTEGRATED – SYSTEM

COMBINE - TREAT - SIMPLIFY



HOW IT WORKS

- ZIP LP PERMITS FREEDOM OF MOVEMENT WITHIN THE SYSTEM AS ARCHWIRES DRIFT POSTERIORLY WITH ANTERIOR DENTAL ALIGNMENT AND OPEN BITE CLOSURE
- NANO HOOKS HAVE BEEN ESPECIALLY DESIGNED ON U2S AND L2S TO ENHANCE SLIDING MECHANICS
- I-ARCH WIRES IS RECOMMENDED FOR EARLIER MOMENTS OF TORQUE TO IMPROVE CONTROL
- TS₂ WORKS AS TONGUE TRAINER AND ANCHORAGE FOR ELASTICS IN PRESENCE OF LATERAL CROSSBITE OR WHEN LATERAL FORCES ARE NECESSARY
- TS₂ AND ZIP LP WORK SYNERGISTICALLY WITH I-ARCH, TO BE CLOSER TO THE CENTER OF RESISTANCE

SIA
ORTHODONTIC
MANUFACTURER

DESIGNING • INNOVATION • MANUFACTURING

info@siaorthodontics.com
www.siaorthodontics.com