

**carefix**<sup>®</sup>  
*Creating functional comfort*



**Passive support  
for personal  
protection**

**Prevents  
hand-to-face  
contamination**

**Barrier mask**  
**Seamless, one-piece protection garment**

**TYTEX**<sup>®</sup>

The Carefix barrier mask is used as a passive support garment that forms a barrier for hand-to-face contamination. It may also help minimize the transmission of airborne organisms from coughing, speaking, sneezing or respiration from persons nearby.

The barrier mask is not a medical supply, but a valid alternative to protect oneself when going outside.

The barrier mask is in a seamless and so-called one-piece tube design, eliminating the use of irritating elastic bands around the ears. The design of the product ensures easy and safe use as well as a good fit, so that the mask stays in place and is extremely comfortable to wear.

Compared to single-use face masks, the Carefix barrier mask is better for the environment, as it is washable and can be reused multiple times.

**Features:**

- Reduces the risk of hand-to-face contamination
- May help minimize the transmission of airborne organisms
- Protects both user and surrounding persons
- Good fit ensures high comfort, correct placement and help avoid too much touching of the mask
- Easy to apply
- Soft, flexible, seamless and breathable material
- Can be folded to small size and easily kept in a bag or pocket
- Washable and reusable - can be washed up to 25 times
- Latex-free material for the safety of the user
- Tested for harmful substances according to Oeko-Tex standard 100

**Recommended for:**

The seamless barrier mask is recommended to be used as a passive support by people who are looking for personal protection as a result of social distancing and other government regulations. The mask is a reusable mask for civil use only.



**Composition:**

94% Polyamide, 6% Elastane

**Washing instructions:**



95°C/203°F

**Available sizes:**

	S/M		M/L	
	cm	in	cm	in
	46-52	18-21	51-60	20-24

Circumference from neck to mid nose.

**Available colours:**



**Packaging:**

- Bags of 5 pcs. per size and colour in boxes of 400 pcs.
- Bags of 100 pcs. per size and colour in boxes of 400 pcs.

**Instructions for use**

Make sure to choose the right size. The mask must sit tight around your mouth and nose.

**Application**

- Wash your hands thoroughly prior to application of the mask.
- Apply the mask from the top of the head with the wide part toward the nose and the narrow part toward the neck.
- With both hands, pull the mask down around the head, making sure to cover the nose and mouth.
- Arrange the mask just under the ears and around the neck. Be careful not to touch mouth, nose or eyes when applying the mask.



**Removal**

- Wash your hands thoroughly before the mask is removed.
- Remove the mask by lifting the narrow part covering the neck over the head and then remove the mask from mouth and nose, while still holding on to the narrow part. Be careful not to touch mouth, nose or eyes when removing the mask.



**IMPORTANT:** Mask is effective only when used in combination with frequent hand-cleaning with alcohol-based hand rub or soap and water. Mask can be reused after washing at 95°C/203°F, maximum 25 times.

Contact Information:  
**JD Honigberg International**  
medical@jdhintl.com  
USA: (1) 847-412-0200



TYTEX A/S has been certified according to ISO 9001, ISO 14001 and ISO 13485

