

03  
**Dr. Classic - Stay Classic Designing**

Fully Automatic  
5 Preset sterilization programs  
Ster. Temp.: 105–134°C  
Ster. Time: 4–60 min  
Dry time adjustable: 1–60min  
B&D / Helix / Vacuum testing program  
Built-in Printer  
Big Stainless Steel Water Tanks  
High quality bacteria filter

**Parameters**

Model	LF-18L	LF-23L
Capacity	18L	23L
Chamber size	Ø2075mm	Ø2075mm
Display	LCD	LCD
Lock	Double Lock	Double Lock
Product size	63 X 50 X 43cm	63 X 50 X 43cm
GW / NetW	61 / 56 KG	63 / 56 KG
Power	220V / 50Hz 2.4 KW	220V / 50Hz 2.4 KW
Output	USB + Printer	Printer

**Function**

Program	121°C Solid	121°C Universal	134°C Solid	134°C Universal	Cotton	BD/Helix Test	Vacuum Test
Vacuum Times	1	3	1	3	3	3	3
Temp (°C)	121	121	134	134	134	134	/
Pressure(KPA)	110	110	210	210	210	210	-80
Time(Min.)	20	20	4	4	4	3.5	15



04  
**LF-18L-E & LF-23L-E**

Fully Automatic  
7 Preset sterilization programs  
1 Custom program  
Ster. Temp.: 105–134°C  
Ster. Time: 4–60 min  
Dry time: 1–25min  
B&D / Helix / Vacuum testing program  
Cleaning Program  
Printer & USB Output  
Big Stainless Steel Water Tanks  
High quality bacteria filter  
User Password

**Parameters**

Model	LF-18L-E	LF-23L-E
Capacity	18L	23L
Chamber size	Ø2075mm	Ø2075mm
Display	LCD	LCD
Lock	Double Lock	Double Lock
Product size	63 X 50 X 43cm	63 X 50 X 43cm
GW / NetW	61 / 56 KG	63 / 57 KG
Power	220V / 50Hz 2.4 KW	220V / 50Hz 2.4 KW
Output	USB + Printer	USB + Printer

**Function**

Program	121°C Solid	121°C Hollow	121°C Porous	121°C Solid	121°C Hollow	121°C Porous	Cotton	Prion	User defined	BD/Helix Test	Vacuum Test	Clean
Vacuum Times	1	3	3	1	3	3	3	3	User defined	3	1	3
Temp (°C)	121	121	121	134	134	134	134	134	User defined	134	/	105
Pressure(KPA)	110	110	110	210	210	210	210	210	User defined	210	-80	20
Time(Min.)	20	20	20	4	4	4	4	19	User defined	3.5	15	5



05  
**Dr. Compac - Fashion Autoclave  
3 - year warranty**

Fast sterilization  
Fully Automatic  
7 Preset Sterilization Program  
Ster. Temp.: 105–134°C  
Ster. Time: 4–60 min  
Dry Time: 1–25 min  
Vacuum Times: 1–10times  
B&D / HELIX / VACUUM Testing  
Cleaning Program  
Printer & USB Output Available  
High quality plastic water tanks  
High-quality Bacteria filter

**Parameters**

Model	LF-18L-G	LF-23L-G
Capacity	18L	23L
Chamber size	Ø2075mm	Ø2075mm
Display	LCD	LCD
Lock	Double Lock	Double Lock
Product size	63 X 50 X 43cm	63 X 50 X 43cm
GW / NetW	61 / 56 KG	63 / 57 KG
Power	220V / 50Hz 2.4 KW	220V / 50Hz 2.4 KW
Output	USB + Printer	USB + Printer

**Function**

Program	121°C Solid	121°C Hollow	121°C Porous	121°C Solid	121°C Hollow	121°C Porous	Cotton	Prion	User defined	BD/Helix Test	Vacuum Test	Clean
Vacuum Times	1	3	3	1	3	3	3	3	User defined	3	1	3
Temp (°C)	121	121	121	134	134	134	134	134	User defined	134	/	105
Pressure(KPA)	110	110	110	210	210	210	210	210	User defined	210	-80	20
Time(Min.)	20	20	20	4	4	4	4	19	User defined	3.5	15	5

